

Product Name:	LF141WWW/0HC
Description:	<p>General Purpose Auto Defrost Hydrocarbon Laboratory Freezer 14 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: -15°-25°C Four adjustable shelves and wire basket Magnetic door gasket for positive seal Heavy duty door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation Hydrocarbon, natural refrigerant (R600a) One year parts and labor warranty; plus additional four year compressor parts warranty Exterior dimensions: 28" W x 26 1/4" D x 67 3/4" H Auto defrost 115V, 60Hz, 5 Amps, 1/4 HP UL/C-UL listed Shipping weight: 185 lbs DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.</p>