



Hydrogen Peroxide, 30% Solution AR (ACS)

Product No. 5240
Specifications current as of: Jun 8
2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Assay (H ₂ O ₂)	29.0 - 32.0 %
Ammonium (NH ₄)	max. 0.0005 %
Chloride (Cl)	max. 0.0003 %
Color (APHA)	max. 10
Heavy Metals (as Pb)	max. 0.0001 %
Iron (Fe)	max. 0.00005 %
Nitrate (NO ₃)	max. 0.0002 %
Phosphate (PO ₄)	max. 0.0002 %
Residue after Evaporation	max. 0.002 %
Sulfate (SO ₄)	max. 0.0005 %
Titration Acid (meq/g)	max. 0.0006

Product Information (not specifications):
Appearance (clear, colorless liquid)

For Laboratory, Research or Manufacturing Use