



317PR

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

A Geno Technology, Inc. (USA) brand name


OneQuant™ HOOK™ NHS- Biotin

(Cat. #786-083)



think proteins! think G-Biosciences www.GBiosciences.com

INTRODUCTION



NHS-Biotin is a commonly used amine reactive biotin-labeling reagent. OneQuant™ HOOK™ NHS-Biotin are single aliquots of NHS-Biotin reagent in convenient, premeasured microtubes. Each microtube contains 2mg of the biotin reagent. Simply dissolve one vial with appropriate solvent and add to the reaction.

Chemical name: N-Hydroxysuccinimidobiotin

Formula: C₁₄H₁₉N₃O₅S

Molecular Weight: 341.38

CAS No: 35013-72-0

Spacer Arm Length: 13.5Å

Reaction pH: 7-9

Membrane permeable

Solvent: DMSO or DMF

ITEMS SUPPLIED

Description	Size
OneQuant™ HOOK™ NHS-Biotin	2mg/vial, 1 strip of 8 vials

STORAGE CONDITIONS

It is shipped at ambient temperature. Upon receipt store product at 4°C protected from moisture. Stable for one year.

RELATED PRODUCTS

Download our Protein Labeling & Conjugation and Antibody Production Handbook.



<http://info2.gbiosciences.com/complete-protein-labeling-conjugation-handbook>

<http://info2.gbiosciences.com/complete-antibody-production-handbook>

For other related products, visit our website at www.GBiosciences.com or contact us.

Last saved: 3/19/2014 TNN



317PR

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

A Geno Technology, Inc. (USA) brand name



www.GBiosciences.com