

## MSDS #8242 COVER SHEET

<b>35070</b>	<b>Pierce® Fast Western Blot, SuperSignal® WEST Dura, Mouse</b>
Component #	Description
1856145	Luminol / Enhancer Solution
1856146	Stable Peroxide Solution
1862119	Antibody Diluent
1862120	10X Wash Buffer
1862115	Mouse Optimized HRP Reagent

<b>35075</b>	<b>Pierce® Fast Western Blot, SuperSignal® WEST Dura, Mouse</b>
Component #	Description
1856157	Luminol / Enhancer Solution
1856158	Stable Peroxide Solution
1862117	Antibody Diluent
1862118	10X Wash Buffer
1862112	Mouse Optimized HRP Reagent

# Material Safety Data Sheet

## 1. Product and company identification

**Product name** : Luminol/Enhancer Solutions

**Synonym** : Supersignal® ELISA Pico Luminol/Enhancer; Supersignal® Stable Luminol/Enhancer Solution; Supersignal® ELISA Femto Luminol/Enhancer; Supersignal® West Femto Luminol/Enhancer Solution; Supersignal® West Pico Luminol/Enhancer Solution; Supersignal® West Dura Luminol/Enhancer Solution; Luminol/Enhancer Solution; North2South Luminol/Enhancer Solution; LightShift® Luminol/Enhancer Solution

**Supplier** : Thermo Fisher Scientific  
Pierce Biotechnology  
P.O. Box 117  
Rockford, IL 61105  
United States  
815.968.0747 or  
800.874.3723

**Manufacturer** : Thermo Fisher Scientific  
Pierce Biotechnology  
P.O. Box 117  
Rockford, IL 61105  
United States  
815.968.0747  
800.874.3723

**Product No.** : 0034075B 0034080B 0034095B 0034095BE 0037070B 0037075B 0089880E 1856136 1856136B 1856145 1856151 1856155 1856157 1856189 1856192 1856193 1856229 1856302 1856304 1856320 1858089 1858114 1858410 1858593 1858783 1859022 1859024 1859027 1859028 1859120 1859675 1859676 1859678 1860852 1860854 1862123 1900252 1900253 1900288 1900289 NCI0025 NCI2012 NCI3408B NCI3409B NCI3708B NCI4095B NCI7075B NCI1858089

**MSDS #** : 7337

**Validation date** : 10/20/2009.

**Print date** : 10/20/2009.

**Responsible name** : **MSDS (Regulatory Affairs)**

**In case of emergency** : CHEMTREC: **Use of** : Refer to the instruction booklet for proper and intended use.  
800.424.9300 **Substance/Preparation** : Otherwise, contact supplier for specific applications.  
OUTSIDE US:  
202.483.7616

## 2. Hazards identification

**Physical state** : Liquid. [Clear sparkling liquid.]

**OSHA/HCS status** : This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).

**Emergency overview** : WARNING !  
CAUSES RESPIRATORY TRACT, EYE AND SKIN IRRITATION. MAY BE HARMFUL IF SWALLOWED. CONTAINS MATERIAL THAT CAN CAUSE TARGET ORGAN DAMAGE.  
Harmful if swallowed. Irritating to eyes, respiratory system and skin. Avoid exposure - obtain special instructions before use. Do not breathe vapor or mist. Do not ingest. Avoid contact with eyes, skin and clothing. Contains material that can cause target organ damage. Use only with adequate ventilation. Keep container tightly closed and sealed until ready for use. Wash thoroughly after handling.

**Routes of entry** : Eye contact. Inhalation. Ingestion.

**Potential acute health effects**

10/20/2009.

Biosciences Division 3747 N. Meridian Rd. Rockford, IL (815) 968-0747 www.thermo.com/pierce  
Pierce P.O. Box 117 61105 (815) 968-7316 Fax

1/9

## Luminol/Enhancer Solutions

## 2. Hazards identification

**Inhalation** : Irritating to respiratory system. Exposure to decomposition products may cause a health hazard. Serious effects may be delayed following exposure.

**Ingestion** : Harmful if swallowed.

**Skin** : Irritating to skin.

**Eyes** : Irritating to eyes.

**Potential chronic health effects**

**Chronic effects** : Contains material that can cause target organ damage.

**Carcinogenicity** : No known significant effects or critical hazards.

**Mutagenicity** : No known significant effects or critical hazards.

**Teratogenicity** : No known significant effects or critical hazards.

**Developmental effects** : No known significant effects or critical hazards.

**Fertility effects** : No known significant effects or critical hazards.

**Target organs** : Contains material which causes damage to the following organs: upper respiratory tract, skin, central nervous system (CNS), eye, lens or cornea.  
Contains material which may cause damage to the following organs: kidneys, the nervous system, liver, cardiovascular system.

**Over-exposure signs/symptoms**

**Inhalation** : Adverse symptoms may include the following:  
respiratory tract irritation  
coughing

**Ingestion** : No specific data.

**Skin** : Adverse symptoms may include the following:  
irritation  
redness

**Eyes** : Adverse symptoms may include the following:  
pain or irritation  
watering  
redness

**Medical conditions aggravated by over-exposure** : Pre-existing digestive disorders and disorders involving any other target organs mentioned in this MSDS as being at risk may be aggravated by over-exposure to this product.

The preparation is not classified as dangerous according to Directive 1999/45/EC and its amendments.

**Classification** : Not classified.

**See toxicological information (section 11)**

## 3. Composition/information on ingredients

### United States

Name	CAS number	%
Tris (Tris(hydroxymethyl)aminomethane)	77-86-1	1 - 3
Ethylene Glycol	107-21-1	1 - 3

### Europe

**Substance/preparation** : Preparation

Ingredient name	CAS number	%	EC number	Classification
-----------------	------------	---	-----------	----------------

10/20/2009.

Biosciences Division 3747 N. Meridian Rd. Rockford, IL (815) 968-0747 www.thermo.com/pierce  
Pierce P.O. Box 117 61105 (815) 968-7316 Fax

2/9

**Lumino/Enhancer Solutions**

**3. Composition/information on ingredients**

Tris (Tris(hydroxymethyl)aminomethane)	77-86-1	1 - 3	201-064-4	Xi; R36/37/38	[1]
Ethylene Glycol	107-21-1	1 - 3	203-473-3	Xn; R22	[1] [2]

There are no ingredients or additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.

[1] Substance classified with a health or environmental hazard

[2] Substance with a workplace exposure limit

Occupational exposure limits, if available, are listed in section 8.

**4. First aid measures**

**Inhalation** : Move exposed person to fresh air. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. Keep person warm and at rest. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation. Get medical attention. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband. In case of inhalation of decomposition products in a fire, symptoms may be delayed. The exposed person may need to be kept under medical surveillance for 48 hours.

**Ingestion** : Wash out mouth with water. Remove dentures if any. Move exposed person to fresh air. Keep person warm and at rest. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Get medical attention. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

**Skin contact** : Flush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Continue to rinse for at least 10 minutes. Get medical attention. Wash clothing before reuse. Clean shoes thoroughly before reuse.

**Eye contact** : Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Get medical attention.

**Protection of first-aiders** : No action shall be taken involving any personal risk or without suitable training. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation.

See section 11 for more detailed information on health effects and symptoms.

**5. Fire-fighting measures**

**Flammability of the product** : In a fire or if heated, a pressure increase will occur and the container may burst.

**Extinguishing media**

**Suitable** : Use an extinguishing agent suitable for the surrounding fire.

**Not suitable** : None known.

**Lumino/Enhancer Solutions**

**5. Fire-fighting measures**

**Special exposure hazards** : Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training.

**Hazardous combustion products** : Decomposition products may include the following materials:  
carbon oxides  
nitrogen oxides

**Special protective equipment for fire-fighters** : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

**6. Accidental release measures**

**Personal precautions** : No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Avoid breathing vapor or mist. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment (see section 8).

**Environmental precautions** : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).

**Large spill** : Stop leak if without risk. Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Wash spillages into an effluent treatment plant or proceed as follows. Contain and collect spillage with non-combustible, absorbent material e.g. sand, earth, vermiculite or diatomaceous earth and place in container for disposal according to local regulations (see section 13). Dispose of via a licensed waste disposal contractor. Contaminated absorbent material may pose the same hazard as the spilled product. Note: see section 1 for emergency contact information and section 13 for waste disposal.

**Small spill** : Stop leak if without risk. Move containers from spill area. Dilute with water and mop up if water-soluble or absorb with an inert dry material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor.

**7. Handling and storage**

**Handling** : Put on appropriate personal protective equipment (see section 8). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Do not get in eyes or on skin or clothing. Do not breathe vapor or mist. Do not ingest. Use only with adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Empty containers retain product residue and can be hazardous. Do not reuse container.

**Storage** : Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see section 10) and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

**Lumino/Enhancer Solutions**

**8 . Exposure controls/personal protection**

**Product name**

**Exposure limits**

**United States**

ethanediol

**ACGIH TLV (United States, 1/2008).**  
C: 100 mg/m<sup>3</sup> Form: Aerosol  
**OSHA PEL 1989 (United States, 3/1989).**  
CEIL: 125 mg/m<sup>3</sup>  
CEIL: 50 ppm

**Europe**

ethanediol

**EU OEL (Europe, 4/2006). Skin Notes: Indicative**  
Short term limit value: 104 mg/m<sup>3</sup> 15 minute(s).  
Short term limit value: 40 ppm 15 minute(s).  
Limit value: 52 mg/m<sup>3</sup> 8 hour(s).  
Limit value: 20 ppm 8 hour(s).

**Consult local authorities for acceptable exposure limits.**

**Recommended monitoring procedures** : If this product contains ingredients with exposure limits, personal, workplace atmosphere or biological monitoring may be required to determine the effectiveness of the ventilation or other control measures and/or the necessity to use respiratory protective equipment.

**Engineering measures** : Use only with adequate ventilation. If user operations generate dust, fumes, gas, vapor or mist, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure to airborne contaminants below any recommended or statutory limits.

**Hygiene measures** : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.

**Personal protection**

**Respiratory** : Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

**Hands** : Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

**Eyes** : Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts.

**Skin** : Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.

**Environmental exposure controls** : Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.

**Lumino/Enhancer Solutions**

**9 . Physical and chemical properties**

**Physical state** : Liquid. [Clear sparkling liquid.]

**Color** : Colorless. Light Pink

**pH** : 9.5

**Solubility** : Easily soluble in the following materials: cold water and hot water.

**10 . Stability and reactivity**

**Chemical stability** : The product is stable. Under normal conditions of storage and use, hazardous polymerization will not occur.

**Conditions to avoid** : No specific data.

**Incompatible materials** : No specific data.

**Hazardous decomposition products** : Under normal conditions of storage and use, hazardous decomposition products should not be produced.

**Possibility of hazardous reactions** : Will not occur.

**11 . Toxicological information**

**United States**

**Acute toxicity**

Product/ingredient name	Result	Species	Dose	Exposure
trometamol	LD50 Intravenous	Rat	1800 mg/kg	-
	LD50 Oral	Rat	5900 mg/kg	-
ethanediol	LD50 Dermal	Rabbit	9530 uL/kg	-
	LD50	Rat	5010 mg/kg	-
	Intraperitoneal			
	LD50 Intravenous	Rat	3260 mg/kg	-
	LD50 Oral	Rat	4700 mg/kg	-
	LD50 Oral	Rat	4000 mg/kg	-
	LD50	Rat	2800 mg/kg	-
	Subcutaneous			
	LD50 Unreported	Rat	13 g/kg	-
	LDLo	Rat	3300 mg/kg	-
Intramuscular				
LDLo	Rat	3300 mg/kg	-	
Intramuscular				
LDLo Intravenous	Rat	2800 mg/kg	-	
TDL0 Oral	Rat	1000 mg/kg	-	
TDL0	Rat	3000 mg/kg	-	
Subcutaneous				
TDL0 Oral	Rat	120 mg/kg	-	
TDL0 Oral	Rat	1110 mg/kg	-	
TDL0 Oral	Rat	5000 mg/kg	-	

**Conclusion/Summary** : Not available.

**Chronic toxicity**

**Conclusion/Summary** : Not available.

**Carcinogenicity**

**Conclusion/Summary** : Not available.

**Classification**

**Luminal/Enhancer Solutions**

**11 . Toxicological information**

Product/ingredient name	ACGIH	IARC	EPA	NIOSH	NTP	OSHA
trometamol	-	-	-	None.	-	None.
ethanediol	A4	-	-	None.	-	None.

**Mutagenicity**

**Conclusion/Summary** : Not available.

**Teratogenicity**

**Conclusion/Summary** : Not available.

**Reproductive toxicity**

**Conclusion/Summary** : Not available.

**Europe**

**Chronic effects** : No known significant effects or critical hazards.

**Carcinogenicity** : No known significant effects or critical hazards.

**Mutagenicity** : No known significant effects or critical hazards.

**Teratogenicity** : No known significant effects or critical hazards.

**Developmental effects** : No known significant effects or critical hazards.

**Fertility effects** : No known significant effects or critical hazards.

**12 . Ecological information**

**Environmental effects** : No known significant effects or critical hazards.

**United States**

**Aquatic ecotoxicity**

**Conclusion/Summary** : Not available.

**Other adverse effects** : No known significant effects or critical hazards.

**13 . Disposal considerations**

**Waste disposal** : The generation of waste should be avoided or minimized wherever possible. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.

**Hazardous waste** : Within the present knowledge of the supplier, this product is not regarded as hazardous waste, as defined by EU Directive 91/689/EEC.

**Disposal should be in accordance with applicable regional, national and local laws and regulations. Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.**

**Luminal/Enhancer Solutions**

**14 . Transport information**

Regulatory information	UN number	Proper shipping name	Classes	PG*
<b>DOT Classification</b>	Not regulated.	-	-	-
<b>IATA-DGR Class</b>	Not available.	Not available.	Not available.	-

PG\* : Packing group

**15 . Regulatory information**

**United States**

**HCS Classification** : Irritating material  
Target organ effects

**U.S. Federal regulations** : **United States inventory (TSCA 8b)**: Not determined.

**SARA 302/304/311/312 extremely hazardous substances**: No products were found.

**SARA 302/304 emergency planning and notification**: No products were found.

**SARA 302/304/311/312 hazardous chemicals**: trometamol; ethanediol

**SARA 311/312 MSDS distribution - chemical inventory - hazard identification**: trometamol: Immediate (acute) health hazard; ethanediol: Immediate (acute) health hazard, Delayed (chronic) health hazard

**Clean Water Act (CWA) 307**: No products were found.

**Clean Water Act (CWA) 311**: No products were found.

**Clean Air Act (CAA) 112 accidental release prevention**: No products were found.

**Clean Air Act (CAA) 112 regulated flammable substances**: No products were found.

**Clean Air Act (CAA) 112 regulated toxic substances**: No products were found.

**SARA 313**

	Product name	CAS number	Concentration
<b>Form R - Reporting requirements</b>	ethanediol	107-21-1	1 - 3
<b>Supplier notification</b>	ethanediol	107-21-1	1 - 3

SARA 313 notifications must not be detached from the MSDS and any copying and redistribution of the MSDS shall include copying and redistribution of the notice attached to copies of the MSDS subsequently redistributed.

**Canada**

**WHMIS (Canada)**

: Class D-2A: Material causing other toxic effects (Very toxic).  
Class D-2B: Material causing other toxic effects (Toxic).

**Canadian lists**

: **CEPA Toxic substances**: None of the components are listed.  
**Canadian ARET**: None of the components are listed.  
**Canadian NPRI**: The following components are listed: Ethylene glycol  
**Alberta Designated Substances**: None of the components are listed.  
**Ontario Designated Substances**: None of the components are listed.  
**Quebec Designated Substances**: None of the components are listed.

**Canada inventory**

: **Canada inventory**: Not determined.

**EU regulations**

**Risk phrases** : This product is not classified according to EU legislation.  
:

Lumino/Enhancer Solutions

## 15 . Regulatory information

### International regulations

**International lists** : Australia inventory (AICS): Not determined.  
 China inventory (IECSC): Not determined.  
 Korea inventory (KECI): Not determined.  
 Philippines inventory (PICCS): Not determined.  
 Japan inventory (ENCS): Not determined.

## 16 . Other information

**Label requirements** : CAUSES RESPIRATORY TRACT, EYE AND SKIN IRRITATION. MAY BE HARMFUL IF SWALLOWED. CONTAINS MATERIAL THAT CAN CAUSE TARGET ORGAN DAMAGE.

**Hazardous Material Information System (U.S.A.)** :

Health	2
Flammability	0
Physical hazards	0

The customer is responsible for determining the PPE code for this material.

**National Fire Protection Association (U.S.A.)** :



**Date of printing** : 10/20/2009.  
**Date of issue** : 10/20/2009.  
**Date of previous issue** : 7/9/2009.  
**Version** : 1.09

Indicates information that has changed from previously issued version.

**Full text of R-phrases referred to in sections 2 and 3 - Europe** : R22- Harmful if swallowed.  
 R36/37/38- Irritating to eyes, respiratory system and skin.

**Full text of classifications referred to in sections 2 and 3 - Europe** : Xn - Harmful  
 Xi - Irritant

### Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

# Material Safety Data Sheet

## 1. Product and company identification

**Product name** : Stable Peroxide Solutions

**Synonym** : SuperSignal® West Dura Stable Peroxide; SuperSignal® West Femto Max Sensitivity Stable Peroxide; SuperSignal® West Pico Stable Peroxide; SuperSignal® ELISA Pico Stable Peroxide; SuperSignal® ELISA Femto Stable Peroxide; Dura Stable Peroxide; North2South Stable Peroxide Solution; LightShift® Stable Peroxide Solution; Stable Peroxide; Chemiluminescent Substrate B

**Supplier** : Thermo Fisher Scientific  
Pierce Biotechnology  
P.O. Box 117  
Rockford, IL 61105  
United States  
815.968.0747 or  
800.874.3723

**Manufacturer** : Thermo Fisher Scientific  
Pierce Biotechnology  
P.O. Box 117  
Rockford, IL 61105  
United States  
815.968.0747  
800.874.3723

**Product No.** : 0034075A 0034095A 0037070A 0037075A 0089880F 1856135 1856135B 1856146 1856150 1856156 1856158 1856190 1856191 1856194 1856230 1856303 1856305 1856321 1858090 1858115 1858409 1858594 1858780 1859023 1859025 1859026 1859029 1859119 1859674 1859677 1859679 1860853 1860855 1862124 1900266 1900290 1900251 NCI0266 NCI2013 NCI2113 NCI3408A NCI3409A NCI3708A NCI4076 NCI4095A NCI7071 NCI7075A NCI1858090

**MSDS #** : 2290

**Validation date** : 10/20/2009.

**Print date** : 10/20/2009.

**Responsible name** : MSDS (Regulatory Affairs)

**In case of emergency** : CHEMTREC: 800.424.9300  
Use of Substance/Preparation OUTSIDE US: 202.483.7616  
Refer to the instruction booklet for proper and intended use. Otherwise, contact supplier for specific applications.

## 2. Hazards identification

**Physical state** : Liquid. [Clear sparkling liquid.]

**OSHA/HCS status** : While this material is not considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200), this MSDS contains valuable information critical to the safe handling and proper use of the product. This MSDS should be retained and available for employees and other users of this product.

**Emergency overview** : NOT EXPECTED TO PRODUCE SIGNIFICANT ADVERSE HEALTH EFFECTS WHEN THE RECOMMENDED INSTRUCTIONS FOR USE ARE FOLLOWED.  
No known significant effects or critical hazards. Avoid prolonged contact with eyes, skin and clothing.

**Routes of entry** : Dermal contact. Eye contact. Inhalation. Ingestion.

**Potential acute health effects**

**Inhalation** : No known significant effects or critical hazards.

**Ingestion** : No known significant effects or critical hazards.

**Skin** : No known significant effects or critical hazards.

10/20/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

1/8

## Stable Peroxide Solutions

## 2. Hazards identification

**Eyes** : No known significant effects or critical hazards.

**Potential chronic health effects**

**Chronic effects** : No known significant effects or critical hazards.

**Carcinogenicity** : No known significant effects or critical hazards.

**Mutagenicity** : No known significant effects or critical hazards.

**Teratogenicity** : No known significant effects or critical hazards.

**Developmental effects** : No known significant effects or critical hazards.

**Fertility effects** : No known significant effects or critical hazards.

**Over-exposure signs/symptoms**

**Inhalation** : No specific data.

**Ingestion** : No specific data.

**Skin** : No specific data.

**Eyes** : No specific data.

**Medical conditions aggravated by over-exposure** : None known.

The preparation is not classified as dangerous according to Directive 1999/45/EC and its amendments.

**Classification** : Not classified.

See toxicological information (section 11)

## 3. Composition/information on ingredients

**Substance/preparation** : Preparation

There are no ingredients or additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.

## 4. First aid measures

**Inhalation** : Move exposed person to fresh air. Keep person warm and at rest. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. Get medical attention if symptoms occur. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

**Ingestion** : Wash out mouth with water. Remove dentures if any. Move exposed person to fresh air. Keep person warm and at rest. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Get medical attention if symptoms occur. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

**Skin contact** : Flush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Get medical attention if symptoms occur. Wash clothing before reuse. Clean shoes thoroughly before reuse.

10/20/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

2/8

**Stable Peroxide Solutions**

**4 . First aid measures**

**Eye contact** : Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Get medical attention if irritation occurs.

**Protection of first-aiders** : No action shall be taken involving any personal risk or without suitable training.

**See section 11 for more detailed information on health effects and symptoms.**

**5 . Fire-fighting measures**

**Flammability of the product** : In a fire or if heated, a pressure increase will occur and the container may burst.

**Extinguishing media**

**Suitable** : Use an extinguishing agent suitable for the surrounding fire.

**Not suitable** : None known.

**Special exposure hazards** : Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training.

**Hazardous combustion products** : No specific data.

**Special protective equipment for fire-fighters** : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

**6 . Accidental release measures**

**Personal precautions** : No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Avoid breathing vapor or mist. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment (see section 8).

**Environmental precautions** : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).

**Large spill** : Stop leak if without risk. Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Wash spillages into an effluent treatment plant or proceed as follows. Contain and collect spillage with non-combustible, absorbent material e.g. sand, earth, vermiculite or diatomaceous earth and place in container for disposal according to local regulations (see section 13). Dispose of via a licensed waste disposal contractor. Contaminated absorbent material may pose the same hazard as the spilled product. Note: see section 1 for emergency contact information and section 13 for waste disposal.

**Small spill** : Stop leak if without risk. Move containers from spill area. Dilute with water and mop up if water-soluble or absorb with an inert dry material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor.

**7 . Handling and storage**

**Handling** : Put on appropriate personal protective equipment (see section 8). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Do not ingest. Avoid contact with eyes, skin and clothing. Avoid breathing vapor or mist. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Empty containers retain product residue and can be hazardous. Do not reuse container.

**Stable Peroxide Solutions**

**7 . Handling and storage**

**Storage** : Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see section 10) and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

**8 . Exposure controls/personal protection**

**Europe**

No exposure limit value known.

**Consult local authorities for acceptable exposure limits.**

**Recommended monitoring procedures** : If this product contains ingredients with exposure limits, personal, workplace atmosphere or biological monitoring may be required to determine the effectiveness of the ventilation or other control measures and/or the necessity to use respiratory protective equipment.

**Engineering measures** : No special ventilation requirements. Good general ventilation should be sufficient to control worker exposure to airborne contaminants. If this product contains ingredients with exposure limits, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure below any recommended or statutory limits.

**Hygiene measures** : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.

**Personal protection**

**Respiratory** : Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

**Hands** : Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

**Eyes** : Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts.

**Skin** : Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.

**Environmental exposure controls** : Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.



**Stable Peroxide Solutions**

**9 . Physical and chemical properties**

- Physical state** : Liquid. [Clear sparkling liquid.]
- Color** : Colorless.
- pH** : 4.9 to 5.1
- Dispersibility properties** : Easily dispersible in the following materials: cold water and hot water.
- Solubility** : Easily soluble in the following materials: cold water and hot water.

**10 . Stability and reactivity**

- Chemical stability** : The product is stable. Under normal conditions of storage and use, hazardous polymerization will not occur.
- Conditions to avoid** : No specific data.
- Incompatible materials** : No specific data.
- Hazardous decomposition products** : Under normal conditions of storage and use, hazardous decomposition products should not be produced.
- Possibility of hazardous reactions** : Will not occur.

**11 . Toxicological information**

United States

Acute toxicity

- Conclusion/Summary** : To the best of our knowledge, the toxicological properties of this product have not been thoroughly investigated.

Chronic toxicity

- Conclusion/Summary** : Not available.

Carcinogenicity

- Conclusion/Summary** : Not available.

Mutagenicity

- Conclusion/Summary** : Not available.

Teratogenicity

- Conclusion/Summary** : Not available.

Reproductive toxicity

- Conclusion/Summary** : Not available.

Europe

- Chronic effects** : No known significant effects or critical hazards.
- Carcinogenicity** : No known significant effects or critical hazards.
- Mutagenicity** : No known significant effects or critical hazards.
- Teratogenicity** : No known significant effects or critical hazards.
- Developmental effects** : No known significant effects or critical hazards.
- Fertility effects** : No known significant effects or critical hazards.

**Stable Peroxide Solutions**

**12 . Ecological information**

- Environmental effects** : No known significant effects or critical hazards.

United States

Aquatic ecotoxicity

- Conclusion/Summary** : Not available.

Other adverse effects

- : No known significant effects or critical hazards.

**13 . Disposal considerations**

- Waste disposal** : The generation of waste should be avoided or minimized wherever possible. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.

- Hazardous waste** : Within the present knowledge of the supplier, this product is not regarded as hazardous waste, as defined by EU Directive 91/689/EEC.

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.

**14 . Transport information**

Regulatory information	UN number	Proper shipping name	Classes	PG*
<b>DOT Classification</b>	Not regulated.	-	-	-
<b>IATA-DGR Class</b>	Not regulated.	-	-	-

PG\* : Packing group

**15 . Regulatory information**

United States

HCS Classification

- : Not regulated.

U.S. Federal regulations

- United States inventory (TSCA 8b)**: All components are listed or exempted.
- SARA 302/304/311/312 extremely hazardous substances**: No products were found.
- SARA 302/304 emergency planning and notification**: No products were found.
- SARA 302/304/311/312 hazardous chemicals**: No products were found.
- SARA 311/312 MSDS distribution - chemical inventory - hazard identification**: No products were found.
- Clean Water Act (CWA) 307**: No products were found.
- Clean Water Act (CWA) 311**: acetic acid
- Clean Air Act (CAA) 112 accidental release prevention**: No products were found.
- Clean Air Act (CAA) 112 regulated flammable substances**: No products were found.
- Clean Air Act (CAA) 112 regulated toxic substances**: No products were found.

Canada

WHMIS (Canada)

- : Not controlled under WHMIS (Canada).

## 15 . Regulatory information

**Canadian lists** : CEPA Toxic substances: None of the components are listed.  
 Canadian ARET: None of the components are listed.  
 Canadian NPRI: None of the components are listed.  
 Alberta Designated Substances: None of the components are listed.  
 Ontario Designated Substances: None of the components are listed.  
 Quebec Designated Substances: None of the components are listed.

**Canada inventory** : **Canada inventory:** All components are listed or exempted.

**EU regulations**

**Risk phrases** : This product is not classified according to EU legislation.

**International regulations**

**International lists** : **Australia inventory (AICS):** All components are listed or exempted.  
**China inventory (IECSC):** All components are listed or exempted.  
**Korea inventory (KECI):** All components are listed or exempted.  
**Philippines inventory (PICCS):** All components are listed or exempted.  
**Japan inventory (ENCS):** All components are listed or exempted.

## 16 . Other information


**Label requirements** : NOT EXPECTED TO PRODUCE SIGNIFICANT ADVERSE HEALTH EFFECTS WHEN THE RECOMMENDED INSTRUCTIONS FOR USE ARE FOLLOWED.

**Hazardous Material Information System (U.S.A.)** :

Health	0
Flammability	0
Physical hazards	0

The customer is responsible for determining the PPE code for this material.

**National Fire Protection Association (U.S.A.)** :



**Date of printing** : 10/20/2009.

**Date of issue** : 10/20/2009.

**Date of previous issue** : 6/4/2009.

**Version** : 1.08

Indicates information that has changed from previously issued version.

[Notice to reader](#)

## 16 . Other information

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

# Material Safety Data Sheet

## 1. Product and company identification

**Product name** : Fast Western Antibody Diluent  
**Supplier** : Thermo Fisher Scientific  
 Pierce Biotechnology  
 P.O. Box 117  
 Rockford, IL 61105  
 United States  
 815.968.0747 or  
 800.874.3723

**Manufacturer** : Thermo Fisher Scientific  
 Pierce Biotechnology  
 P.O. Box 117  
 Rockford, IL 61105  
 United States  
 815.968.0747  
 800.874.3723

**Product No.** : 1861766 1861901 1862117 1862119 1862121  
**MSDS #** : 8162  
**Validation date** : 4/10/2009.  
**Print date** : 4/10/2009.  
**Responsible name** : MSDS (Regulatory Affairs)  
**In case of emergency** : CHEMTREC: 800.424.9300  
 OUTSIDE US: 202.483.7616

**Use of Substance/Preparation** : Refer to the instruction booklet for proper and intended use. Otherwise, contact supplier for specific applications.

## 2. Hazards identification

**Physical state** : Liquid.  
**Odor** : Characteristic.  
**OSHA/HCS status** : This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).  
**Emergency overview** : CONTAINS MATERIAL THAT CAN CAUSE TARGET ORGAN DAMAGE.  
 Avoid exposure - obtain special instructions before use. Do not breathe vapor or mist. Contains material that can cause target organ damage.

**Potential acute health effects**

**Inhalation** : No known significant effects or critical hazards.  
**Ingestion** : No known significant effects or critical hazards.  
**Skin** : No known significant effects or critical hazards.  
**Eyes** : No known significant effects or critical hazards.

**Potential chronic health effects**

**Chronic effects** : Contains material that can cause target organ damage.  
**Carcinogenicity** : No known significant effects or critical hazards.  
**Mutagenicity** : No known significant effects or critical hazards.  
**Teratogenicity** : No known significant effects or critical hazards.  
**Developmental effects** : No known significant effects or critical hazards.  
**Fertility effects** : No known significant effects or critical hazards.

4/10/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
P.O. Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

1/8

### Fast Western Antibody Diluent

## 2. Hazards identification

**Target organs** : Contains material which causes damage to the following organs: upper respiratory tract, eye, lens or cornea.  
 Contains material which may cause damage to the following organs: teeth.

**Over-exposure signs/symptoms**

**Inhalation** : No specific data.  
**Ingestion** : No specific data.  
**Skin** : No specific data.  
**Eyes** : No specific data.  
**Medical conditions aggravated by over-exposure** : Pre-existing disorders involving any target organs mentioned in this MSDS as being at risk may be aggravated by over-exposure to this product.

The preparation is not classified as dangerous according to Directive 1999/45/EC and its amendments.

**Classification** : Not classified.

See toxicological information (section 11)

## 3. Composition/information on ingredients

### United States

Name	CAS number	%
Sucrose, 99+% A.C.S. Reagent	57-50-1	3 - 5

### Europe

**Substance/preparation** : Preparation

Ingredient name	CAS number	%	EC number	Classification
Sucrose, 99+% A.C.S. Reagent	57-50-1	3 - 5	200-334-9	Not classified. [2]

There are no ingredients or additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.

[1] Substance classified with a health or environmental hazard

[2] Substance with a workplace exposure limit

Occupational exposure limits, if available, are listed in section 8.

## 4. First aid measures

**Inhalation** : Move exposed person to fresh air. Keep person warm and at rest. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation. Get medical attention. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

4/10/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
P.O. Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

2/8

**Fast Western Antibody Diluent**

**4. First aid measures**

- Ingestion** : Wash out mouth with water. Remove dentures if any. Move exposed person to fresh air. Keep person warm and at rest. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Get medical attention. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.
- Skin contact** : Flush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Continue to rinse for at least 10 minutes. Get medical attention if symptoms occur. Wash clothing before reuse. Clean shoes thoroughly before reuse.
- Eye contact** : Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Get medical attention if irritation occurs.
- Protection of first-aiders** : No action shall be taken involving any personal risk or without suitable training. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation.

See section 11 for more detailed information on health effects and symptoms.

**5. Fire-fighting measures**

- Flammability of the product** : In a fire or if heated, a pressure increase will occur and the container may burst.
- Extinguishing media**
- Suitable** : Use an extinguishing agent suitable for the surrounding fire.
- Not suitable** : None known.
- Special exposure hazards** : Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training.
- Hazardous combustion products** : Decomposition products may include the following materials: carbon oxides
- Special protective equipment for fire-fighters** : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

**6. Accidental release measures**

- Personal precautions** : No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Avoid breathing vapor or mist. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment (see section 8).
- Environmental precautions** : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).
- Large spill** : Stop leak if without risk. Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Wash spillages into an effluent treatment plant or proceed as follows. Contain and collect spillage with non-combustible, absorbent material e.g. sand, earth, vermiculite or diatomaceous earth and place in container for disposal according to local regulations (see section 13). Dispose of via a licensed waste disposal contractor. Contaminated absorbent material may pose the same hazard as the spilled product. Note: see section 1 for emergency contact information and section 13 for waste disposal.

4/10/2009.

3/8

**Fast Western Antibody Diluent**

**6. Accidental release measures**

- Small spill** : Stop leak if without risk. Move containers from spill area. Dilute with water and mop up if water-soluble or absorb with an inert dry material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor.

**7. Handling and storage**

- Handling** : Put on appropriate personal protective equipment (see section 8). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Do not get in eyes or on skin or clothing. Do not breathe vapor or mist. Do not ingest. If during normal use the material presents a respiratory hazard, use only with adequate ventilation or wear appropriate respirator. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Empty containers retain product residue and can be hazardous. Do not reuse container.
- Storage** : Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see section 10) and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

**8. Exposure controls/personal protection**

<u>Product name</u>	<u>Exposure limits</u>
<u>United States</u> sucrose	<b>ACGIH (United States).</b> TWA: 10 mg/m <sup>3</sup> <b>NIOSH (United States).</b> TWA: 5 mg/m <sup>3</sup> <b>ACGIH TLV (United States, 1/2008).</b> TWA: 10 mg/m <sup>3</sup> 8 hour(s). <b>NIOSH REL (United States, 12/2001).</b> TWA: 5 mg/m <sup>3</sup> 10 hour(s). Form: Respirable fraction TWA: 10 mg/m <sup>3</sup> 10 hour(s). Form: Total dust <b>OSHA PEL (United States, 11/2006).</b> TWA: 5 mg/m <sup>3</sup> 8 hour(s). Form: Respirable fraction TWA: 15 mg/m <sup>3</sup> 8 hour(s). Form: Total dust <b>OSHA PEL 1989 (United States, 3/1989).</b> TWA: 5 mg/m <sup>3</sup> 8 hour(s). Form: Respirable fraction TWA: 15 mg/m <sup>3</sup> 8 hour(s). Form: Total dust <b>OSHA PEL (United States). Notes: Respirable</b> TWA: 15 mg/m <sup>3</sup> 8 hour(s).
<u>Europe</u> sucrose	<b>ACGIH (United States).</b> TWA: 10 mg/m <sup>3</sup> <b>ACGIH TLV (United States, 1/2008).</b> TWA: 10 mg/m <sup>3</sup> 8 hour(s).

Consult local authorities for acceptable exposure limits.

- Recommended monitoring procedures** : If this product contains ingredients with exposure limits, personal, workplace atmosphere or biological monitoring may be required to determine the effectiveness of the ventilation or other control measures and/or the necessity to use respiratory protective equipment.

4/10/2009.

4/8

**Fast Western Antibody Diluent**

**8 . Exposure controls/personal protection**

- Engineering measures** : If user operations generate dust, fumes, gas, vapor or mist, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure to airborne contaminants below any recommended or statutory limits.
- Hygiene measures** : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.
- Personal protection**
- Respiratory** : Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.
- Hands** : Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.
- Eyes** : Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts.
- Skin** : Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.
- Environmental exposure controls** : Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.

**9 . Physical and chemical properties**

- Physical state** : Liquid.
- Color** : Amber.
- Odor** : Characteristic.
- pH** : 6 to 8
- Solubility** : Soluble in the following materials: cold water.

**10 . Stability and reactivity**

- Chemical stability** : The product is stable. Under normal conditions of storage and use, hazardous polymerization will not occur.
- Conditions to avoid** : No specific data.
- Incompatible materials** : No specific data.
- Hazardous decomposition products** : Under normal conditions of storage and use, hazardous decomposition products should not be produced.

**Fast Western Antibody Diluent**

**11 . Toxicological information**

**United States**

- Acute toxicity**
- | Product/ingredient name | Result    | Species | Dose        | Exposure |
|-------------------------|-----------|---------|-------------|----------|
| sucrose                 | LD50 Oral | Rat     | 29700 mg/kg | -        |
- Conclusion/Summary** : Not available.
- Chronic toxicity**
- Conclusion/Summary** : Not available.
- Carcinogenicity**
- Conclusion/Summary** : Not available.
- Classification**
- | Product/ingredient name | ACGIH | IARC | EPA | NIOSH | NTP | OSHA |
|-------------------------|-------|------|-----|-------|-----|------|
| sucrose                 | A4    | -    | -   | None. | -   | -    |
- Mutagenicity**
- Conclusion/Summary** : Not available.
- Teratogenicity**
- Conclusion/Summary** : Not available.
- Reproductive toxicity**
- Conclusion/Summary** : Not available.

**Europe**

- Chronic effects** : No known significant effects or critical hazards.
- Carcinogenicity** : No known significant effects or critical hazards.
- Mutagenicity** : No known significant effects or critical hazards.
- Teratogenicity** : No known significant effects or critical hazards.
- Developmental effects** : No known significant effects or critical hazards.
- Fertility effects** : No known significant effects or critical hazards.

**12 . Ecological information**

- Environmental effects** : No known significant effects or critical hazards.
- United States**
- Aquatic ecotoxicity**
- Conclusion/Summary** : Not available.

**13 . Disposal considerations**

- Waste disposal** : The generation of waste should be avoided or minimized wherever possible. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.
- Hazardous waste** : Within the present knowledge of the supplier, this product is not regarded as hazardous waste, as defined by EU Directive 91/689/EEC.
- Disposal should be in accordance with applicable regional, national and local laws and regulations.**

Fast Western Antibody Diluent

### 13 . Disposal considerations

Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.

### 14 . Transport information

Regulatory information	UN number	Proper shipping name	Classes	PG*
DOT Classification	Not regulated.	-	-	-
IATA-DGR Class	Not regulated.	-	-	-

PG\* : Packing group

### 15 . Regulatory information

#### United States

- HCS Classification** : Target organ effects  
**TSCA 8(a) IUR:** sucrose  
**U.S. Federal regulations** : TSCA 8(a) IUR: sucrose  
**United States inventory (TSCA 8b):** Not determined.  
**SARA 302/304/311/312 extremely hazardous substances:** No products were found.  
**SARA 302/304 emergency planning and notification:** No products were found.  
**SARA 302/304/311/312 hazardous chemicals:** sucrose  
**SARA 311/312 MSDS distribution - chemical inventory - hazard identification:**  
 sucrose: Delayed (chronic) health hazard  
**Clean Water Act (CWA) 307:** No products were found.  
**Clean Water Act (CWA) 311:** hydrochloric acid  
**Clean Air Act (CAA) 112 accidental release prevention:** hydrochloric acid  
**Clean Air Act (CAA) 112 regulated flammable substances:** hydrochloric acid  
**Clean Air Act (CAA) 112 regulated toxic substances:** hydrochloric acid

#### Canada

- WHMIS (Canada)** : Not controlled under WHMIS (Canada).  
**Canadian lists** : **CEPA Toxic substances:** None of the components are listed.  
**Canadian ARET:** None of the components are listed.  
**Canadian NPRI:** None of the components are listed.  
**Alberta Designated Substances:** None of the components are listed.  
**Ontario Designated Substances:** None of the components are listed.  
**Quebec Designated Substances:** None of the components are listed.  
**Canada inventory:** Not determined.

#### EU regulations

- Risk phrases** : This product is not classified according to EU legislation.  
 :

#### International regulations

4/10/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

7/8

Fast Western Antibody Diluent

### 15 . Regulatory information

- International lists** : **Australia inventory (AICS):** Not determined.  
**China inventory (IECSC):** Not determined.  
**Korea inventory (KECI):** Not determined.  
**Philippines inventory (PICCS):** Not determined.  
**Japan inventory (ENCS):** Not determined.

### 16 . Other information


- Label requirements** : CONTAINS MATERIAL THAT CAN CAUSE TARGET ORGAN DAMAGE.

**Hazardous Material Information System (U.S.A.)** :

Health	1
Flammability	0
Physical hazards	0

The customer is responsible for determining the PPE code for this material.

**National Fire Protection Association (U.S.A.)** :



**Date of printing** : 4/10/2009.

**Date of issue** : 4/10/2009.

**Date of previous issue** : 1/23/2009.

**Version** : 1.01

Indicates information that has changed from previously issued version.

#### Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

4/10/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

8/8

# Material Safety Data Sheet

## 1. Product and company identification

**Product name** : Fast Western 10X Wash Buffer

**Supplier** : Thermo Fisher Scientific  
Pierce Biotechnology  
P.O. Box 117  
Rockford, IL 61105  
United States  
815.968.0747 or  
800.874.3723

**Manufacturer** : Thermo Fisher Scientific  
Pierce Biotechnology  
P.O. Box 117  
Rockford, IL 61105  
United States  
815.968.0747  
800.874.3723

**Product No.** : 1861767 1861902 1862118 1862120 1862122

**MSDS #** : 8163

**Validation date** : 4/10/2009.

**Print date** : 4/10/2009.

**Responsible name** : **MSDS (Regulatory Affairs)**

**In case of emergency** : CHEMTREC: **Use of Substance/Preparation**  
800.424.9300  
OUTSIDE US:  
202.483.7616

Refer to the instruction booklet for proper and intended use. Otherwise, contact supplier for specific applications.

## 2. Hazards identification

**Physical state** : Liquid.

**Odor** : Characteristic.

**OSHA/HCS status** : This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).

**Emergency overview** : WARNING !  
CAUSES RESPIRATORY TRACT, EYE AND SKIN IRRITATION. CONTAINS MATERIAL THAT CAN CAUSE TARGET ORGAN DAMAGE.  
Irritating to eyes, respiratory system and skin. Avoid exposure - obtain special instructions before use. Do not breathe vapor or mist. Avoid contact with eyes, skin and clothing. Contains material that can cause target organ damage. Use only with adequate ventilation. Keep container tightly closed and sealed until ready for use. Wash thoroughly after handling.

**Routes of entry** : Dermal contact. Eye contact. Inhalation. Ingestion.

**Potential acute health effects**

**Inhalation** : Irritating to respiratory system. Exposure to decomposition products may cause a health hazard. Serious effects may be delayed following exposure.

**Ingestion** : No known significant effects or critical hazards.

**Skin** : Irritating to skin.

**Eyes** : Irritating to eyes.

**Potential chronic health effects**

**Chronic effects** : Contains material that can cause target organ damage.

**Carcinogenicity** : No known significant effects or critical hazards.

4/10/2009.

1/8

Fast Western 10X Wash Buffer

## 2. Hazards identification

**Mutagenicity** : No known significant effects or critical hazards.

**Teratogenicity** : No known significant effects or critical hazards.

**Developmental effects** : No known significant effects or critical hazards.

**Fertility effects** : No known significant effects or critical hazards.

**Target organs** : Contains material which causes damage to the following organs: kidneys, liver. Contains material which may cause damage to the following organs: skin, eyes, stomach.

**Over-exposure signs/symptoms**

**Inhalation** : Adverse symptoms may include the following:  
respiratory tract irritation  
coughing

**Ingestion** : No specific data.

**Skin** : Adverse symptoms may include the following:  
irritation  
redness

**Eyes** : Adverse symptoms may include the following:  
pain or irritation  
watering  
redness

**Medical conditions aggravated by over-exposure** : Pre-existing disorders involving any target organs mentioned in this MSDS as being at risk may be aggravated by over-exposure to this product.

The preparation is not classified as dangerous according to Directive 1999/45/EC and its amendments.

**Classification** : Not classified.

See toxicological information (section 11)

## 3. Composition/information on ingredients

### United States

Name	CAS number	%
Sodium Chloride	7647-14-5	7 - 10
Tris (Tris(hydroxymethyl)aminomethane)	77-86-1	1 - 3

**Substance/preparation** : Preparation

**There are no ingredients or additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.**

## 4. First aid measures

**Inhalation** : Move exposed person to fresh air. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. Keep person warm and at rest. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation. Get medical attention. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband. In case of inhalation of decomposition products in a fire, symptoms may be delayed.

4/10/2009.

2/8

#### 4. First aid measures

The exposed person may need to be kept under medical surveillance for 48 hours.

- Ingestion** : Wash out mouth with water. Remove dentures if any. Move exposed person to fresh air. Keep person warm and at rest. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Get medical attention. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.
- Skin contact** : Flush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Continue to rinse for at least 10 minutes. Get medical attention. Wash clothing before reuse. Clean shoes thoroughly before reuse.
- Eye contact** : Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Get medical attention.
- Protection of first-aiders** : No action shall be taken involving any personal risk or without suitable training. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation.

See section 11 for more detailed information on health effects and symptoms.

#### 5. Fire-fighting measures

- Flammability of the product** : In a fire or if heated, a pressure increase will occur and the container may burst.
- Extinguishing media**
- Suitable** : Use an extinguishing agent suitable for the surrounding fire.
- Not suitable** : None known.
- Special exposure hazards** : Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training.
- Hazardous combustion products** : Decomposition products may include the following materials:  
carbon oxides  
nitrogen oxides  
halogenated compounds  
metal oxide/oxides
- Special protective equipment for fire-fighters** : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

#### 6. Accidental release measures

- Personal precautions** : No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Avoid breathing vapor or mist. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment (see section 8).
- Environmental precautions** : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).

#### 6. Accidental release measures

- Large spill** : Stop leak if without risk. Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Wash spillages into an effluent treatment plant or proceed as follows. Contain and collect spillage with non-combustible, absorbent material e.g. sand, earth, vermiculite or diatomaceous earth and place in container for disposal according to local regulations (see section 13). Dispose of via a licensed waste disposal contractor. Contaminated absorbent material may pose the same hazard as the spilled product. Note: see section 1 for emergency contact information and section 13 for waste disposal.
- Small spill** : Stop leak if without risk. Move containers from spill area. Dilute with water and mop up if water-soluble or absorb with an inert dry material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor.

#### 7. Handling and storage

- Handling** : Put on appropriate personal protective equipment (see section 8). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Do not get in eyes or on skin or clothing. Do not breathe vapor or mist. Do not ingest. Use only with adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Empty containers retain product residue and can be hazardous. Do not reuse container.
- Storage** : Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see section 10) and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

#### 8. Exposure controls/personal protection

- Europe**  
No exposure limit value known.
- Consult local authorities for acceptable exposure limits.**
- Recommended monitoring procedures** : If this product contains ingredients with exposure limits, personal, workplace atmosphere or biological monitoring may be required to determine the effectiveness of the ventilation or other control measures and/or the necessity to use respiratory protective equipment.
- Engineering measures** : Use only with adequate ventilation. If user operations generate dust, fumes, gas, vapor or mist, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure to airborne contaminants below any recommended or statutory limits.
- Hygiene measures** : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.
- Personal protection**
- Respiratory** : Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.



**Fast Western 10X Wash Buffer**

**8 . Exposure controls/personal protection**

- Hands** : Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.
- Eyes** : Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts.
- Skin** : Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.
- Environmental exposure controls** : Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.

**9 . Physical and chemical properties**

- Physical state** : Liquid.
- Color** : Colorless.
- Odor** : Characteristic.
- pH** : 6 to 8
- Solubility** : Soluble in the following materials: cold water.

**10 . Stability and reactivity**

- Chemical stability** : The product is stable. Under normal conditions of storage and use, hazardous polymerization will not occur.
- Conditions to avoid** : No specific data.
- Incompatible materials** : No specific data.
- Hazardous decomposition products** : Under normal conditions of storage and use, hazardous decomposition products should not be produced.

**11 . Toxicological information**

United States

Acute toxicity

Product/ingredient name	Result	Species	Dose	Exposure
sodium chloride	LD50 Dermal	Rabbit	>10 gm/kg	-
	LD50 Dermal	Rabbit	10000 mg/kg	-
	LD50	Rat	2600 mg/kg	-
	Intraperitoneal			
LD50 Oral	Rat	3000 mg/kg	-	
	LDLo	Rat	3.72 gm/kg	-
Intraperitoneal	Rat			
	LDLo	Rat	3500 mg/kg	-
Subcutaneous				
LC50 Inhalation	Rat	>42 gm/m3	1 hours	
Dusts and mists				
1,3-propanediol, 2-amino-2-	LD50 Intravenous	Rat	1800 mg/kg	-

4/10/2009.

5/8

**Fast Western 10X Wash Buffer**

**11 . Toxicological information**

(hydroxymethyl)- LD50 Oral Rat 5900 mg/kg -

**Conclusion/Summary** : Not available.

Chronic toxicity

**Conclusion/Summary** : Not available.

Carcinogenicity

**Conclusion/Summary** : Not available.

Classification

Product/ingredient name	ACGIH	IARC	EPA	NIOSH	NTP	OSHA
sodium chloride	-	-	-	None.	-	None.
1,3-propanediol, 2-amino-2-(hydroxymethyl)-	-	-	-	None.	-	None.

Mutagenicity

**Conclusion/Summary** : Not available.

Teratogenicity

**Conclusion/Summary** : Not available.

Reproductive toxicity

**Conclusion/Summary** : Not available.

Europe

- Chronic effects** : No known significant effects or critical hazards.
- Carcinogenicity** : No known significant effects or critical hazards.
- Mutagenicity** : No known significant effects or critical hazards.
- Teratogenicity** : No known significant effects or critical hazards.
- Developmental effects** : No known significant effects or critical hazards.
- Fertility effects** : No known significant effects or critical hazards.

**12 . Ecological information**

**Environmental effects** : No known significant effects or critical hazards.

United States

Aquatic ecotoxicity

**Conclusion/Summary** : Not available.

**13 . Disposal considerations**

**Waste disposal** : The generation of waste should be avoided or minimized wherever possible. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.

**Hazardous waste** : Within the present knowledge of the supplier, this product is not regarded as hazardous waste, as defined by EU Directive 91/689/EEC.

**Disposal should be in accordance with applicable regional, national and local laws and regulations.**

4/10/2009.

6/8

Fast Western 10X Wash Buffer

### 13 . Disposal considerations

Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.

### 14 . Transport information

Regulatory information	UN number	Proper shipping name	Classes	PG*
DOT Classification	Not regulated.	-	-	-
IATA-DGR Class	Not regulated.	-	-	-

PG\* : Packing group

### 15 . Regulatory information

#### United States

**HCS Classification** : Irritating material  
Target organ effects

**U.S. Federal regulations** : TSCA 8(a) PAIR: polyethylene glycol octaphenyl ether  
**United States inventory (TSCA 8b)**: All components are listed or exempted.  
**SARA 302/304/311/312 extremely hazardous substances**: No products were found.  
**SARA 302/304 emergency planning and notification**: No products were found.  
**SARA 302/304/311/312 hazardous chemicals**: sodium chloride  
**SARA 311/312 MSDS distribution - chemical inventory - hazard identification**: sodium chloride: Immediate (acute) health hazard, Delayed (chronic) health hazard  
**Clean Water Act (CWA) 307**: No products were found.  
**Clean Water Act (CWA) 311**: No products were found.  
**Clean Air Act (CAA) 112 accidental release prevention**: ethylene oxide  
**Clean Air Act (CAA) 112 regulated flammable substances**: No products were found.  
**Clean Air Act (CAA) 112 regulated toxic substances**: ethylene oxide

#### California Prop. 65

**WARNING:** This product contains a chemical known to the State of California to cause cancer and birth defects or other reproductive harm.

Ingredient name	Cancer	Reproductive	No significant risk level	Maximum acceptable dosage level
ethylene oxide	Yes.	Yes.	Yes.	Yes.

#### Canada

**WHMIS (Canada)** : Class D-2B: Material causing other toxic effects (Toxic).  
**Canada inventory** : **Canada inventory**: At least one component is not listed.

#### EU regulations

**Risk phrases** : This product is not classified according to EU legislation.  
:

#### International regulations

4/10/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

7/8

Fast Western 10X Wash Buffer

### 15 . Regulatory information

**International lists** : **Australia inventory (AICS)**: All components are listed or exempted.  
**China inventory (IECSC)**: All components are listed or exempted.  
**Korea inventory (KECI)**: All components are listed or exempted.  
**Philippines inventory (PICCS)**: All components are listed or exempted.  
**Japan inventory (ENCS)**: All components are listed or exempted.

### 16 . Other information


**Label requirements** : CAUSES RESPIRATORY TRACT, EYE AND SKIN IRRITATION. CONTAINS MATERIAL THAT CAN CAUSE TARGET ORGAN DAMAGE.

**Hazardous Material Information System (U.S.A.)** :

Health	2
Flammability	0
Physical hazards	0

The customer is responsible for determining the PPE code for this material.

**National Fire Protection Association (U.S.A.)** :



**Date of printing** : 4/10/2009.  
**Date of issue** : 4/10/2009.  
**Date of previous issue** : 3/17/2009.  
**Version** : 1.02

Indicates information that has changed from previously issued version.

#### Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

4/10/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

8/8

# Material Safety Data Sheet

## 1. Product and company identification

<b>Product name</b>	: Mouse Optimized HRP Reagent, Dura		
<b>Supplier</b>	: Thermo Fisher Scientific Pierce Biotechnology P.O. Box 117 Rockford, IL 61105 United States 815.968.0747 or 800.874.3723	<b>Manufacturer</b>	: Thermo Fisher Scientific Pierce Biotechnology P.O. Box 117 Rockford, IL 61105 United States 815.968.0747 800.874.3723
<b>Product No.</b>	: 1862112 1862115 1901768		
<b>MSDS #</b>	: 8243		
<b>Validation date</b>	: 10/29/2009.		
<b>Print date</b>	: 10/29/2009.		
<b>Responsible name</b>	: <b>MSDS Specialist</b>		
<b>In case of emergency</b>	: CHEMTREC: 800.424.9300 OUTSIDE US: 202.483.7616	<b>Use of Substance/Preparation</b>	Refer to the instruction booklet for proper and intended use. Otherwise, contact supplier for specific applications.

## 2. Hazards identification

<b>Physical state</b>	: Liquid.
<b>Odor</b>	: Characteristic.
<b>OSHA/HCS status</b>	: While this material is not considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200), this MSDS contains valuable information critical to the safe handling and proper use of the product. This MSDS should be retained and available for employees and other users of this product.
<b>Emergency overview</b>	: NOT EXPECTED TO PRODUCE SIGNIFICANT ADVERSE HEALTH EFFECTS WHEN THE RECOMMENDED INSTRUCTIONS FOR USE ARE FOLLOWED.  No known significant effects or critical hazards. Avoid prolonged contact with eyes, skin and clothing.
<b>Routes of entry</b>	: Dermal contact. Eye contact. Inhalation. Ingestion.
<b>Potential acute health effects</b>	
<b>Inhalation</b>	: No known significant effects or critical hazards.
<b>Ingestion</b>	: No known significant effects or critical hazards.
<b>Skin</b>	: No known significant effects or critical hazards.
<b>Eyes</b>	: No known significant effects or critical hazards.
<b>Potential chronic health effects</b>	
<b>Chronic effects</b>	: No known significant effects or critical hazards.
<b>Carcinogenicity</b>	: No known significant effects or critical hazards.
<b>Mutagenicity</b>	: No known significant effects or critical hazards.
<b>Teratogenicity</b>	: No known significant effects or critical hazards.
<b>Developmental effects</b>	: No known significant effects or critical hazards.

10/29/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

1/7

**Mouse Optimized HRP Reagent, Dura**

## 2. Hazards identification

<b>Fertility effects</b>	: No known significant effects or critical hazards.
<b>Over-exposure signs/symptoms</b>	
<b>Inhalation</b>	: No specific data.
<b>Ingestion</b>	: No specific data.
<b>Skin</b>	: No specific data.
<b>Eyes</b>	: No specific data.
<b>Medical conditions aggravated by over-exposure</b>	: None known.

The preparation is not classified as dangerous according to Directive 1999/45/EC and its amendments.

**Classification** : Not classified.

**See toxicological information (section 11)**

## 3. Composition/information on ingredients

**Substance/preparation** : Preparation

**There are no ingredients or additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.**

## 4. First aid measures

<b>Inhalation</b>	: Move exposed person to fresh air. Keep person warm and at rest. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. Get medical attention if symptoms occur. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.
<b>Ingestion</b>	: Wash out mouth with water. Remove dentures if any. Move exposed person to fresh air. Keep person warm and at rest. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Get medical attention if symptoms occur. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.
<b>Skin contact</b>	: Flush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Get medical attention if symptoms occur. Wash clothing before reuse. Clean shoes thoroughly before reuse.
<b>Eye contact</b>	: Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Get medical attention if irritation occurs.
<b>Protection of first-aiders</b>	: No action shall be taken involving any personal risk or without suitable training.

**See section 11 for more detailed information on health effects and symptoms.**

10/29/2009.

Biosciences Division  
Pierce

3747 N. Meridian Rd.  
PO Box 117

Rockford, IL  
61105

(815) 968-0747  
(815) 968-7316 Fax

www.thermo.com/pierce

2/7

Mouse Optimized HRP Reagent, Dura

## 5. Fire-fighting measures

**Flammability of the product** : In a fire or if heated, a pressure increase will occur and the container may burst.

### Extinguishing media

**Suitable** : Use an extinguishing agent suitable for the surrounding fire.

**Not suitable** : None known.

**Special exposure hazards** : Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training.

**Hazardous combustion products** : No specific data.

**Special protective equipment for fire-fighters** : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

## 6. Accidental release measures

**Personal precautions** : No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Avoid breathing vapor or mist. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment (see section 8).

**Environmental precautions** : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).

**Large spill** : Stop leak if without risk. Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Wash spillages into an effluent treatment plant or proceed as follows. Contain and collect spillage with non-combustible, absorbent material e.g. sand, earth, vermiculite or diatomaceous earth and place in container for disposal according to local regulations (see section 13). Dispose of via a licensed waste disposal contractor. Contaminated absorbent material may pose the same hazard as the spilled product. Note: see section 1 for emergency contact information and section 13 for waste disposal.

**Small spill** : Stop leak if without risk. Move containers from spill area. Dilute with water and mop up if water-soluble or absorb with an inert dry material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor.

## 7. Handling and storage

**Handling** : Put on appropriate personal protective equipment (see section 8). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Do not ingest. Avoid contact with eyes, skin and clothing. Avoid breathing vapor or mist. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Empty containers retain product residue and can be hazardous. Do not reuse container.

**Storage** : Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see section 10) and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

Mouse Optimized HRP Reagent, Dura

## 8. Exposure controls/personal protection

### Europe

No exposure limit value known.

**Consult local authorities for acceptable exposure limits.**

**Recommended monitoring procedures** : If this product contains ingredients with exposure limits, personal, workplace atmosphere or biological monitoring may be required to determine the effectiveness of the ventilation or other control measures and/or the necessity to use respiratory protective equipment.

**Engineering measures** : No special ventilation requirements. Good general ventilation should be sufficient to control worker exposure to airborne contaminants. If this product contains ingredients with exposure limits, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure below any recommended or statutory limits.

**Hygiene measures** : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.

### Personal protection

**Respiratory** : Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

**Hands** : Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

**Eyes** : Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts.

**Skin** : Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.

**Environmental exposure controls** : Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.

## 9. Physical and chemical properties

**Physical state** : Liquid.

**Color** : Amber.

**Odor** : Characteristic.

**pH** : 6 to 8

**Solubility** : Soluble in the following materials: cold water.

**Mouse Optimized HRP Reagent, Dura**

## 10 . Stability and reactivity

- Chemical stability** : The product is stable. Under normal conditions of storage and use, hazardous polymerization will not occur.
- Conditions to avoid** : No specific data.
- Incompatible materials** : No specific data.
- Hazardous decomposition products** : Under normal conditions of storage and use, hazardous decomposition products should not be produced.

## 11 . Toxicological information

### United States

#### Acute toxicity

- Conclusion/Summary** : To the best of our knowledge, the toxicological properties of this product have not been thoroughly investigated.

#### Chronic toxicity

- Conclusion/Summary** : Not available.

#### Carcinogenicity

- Conclusion/Summary** : Not available.

#### Mutagenicity

- Conclusion/Summary** : Not available.

#### Teratogenicity

- Conclusion/Summary** : Not available.

#### Reproductive toxicity

- Conclusion/Summary** : Not available.

### Europe

- Chronic effects** : No known significant effects or critical hazards.
- Carcinogenicity** : No known significant effects or critical hazards.
- Mutagenicity** : No known significant effects or critical hazards.
- Teratogenicity** : No known significant effects or critical hazards.
- Developmental effects** : No known significant effects or critical hazards.
- Fertility effects** : No known significant effects or critical hazards.

## 12 . Ecological information

- Environmental effects** : No known significant effects or critical hazards.

### United States

#### Aquatic ecotoxicity

- Conclusion/Summary** : Not available.

#### Other adverse effects

- Conclusion/Summary** : No known significant effects or critical hazards.

**Mouse Optimized HRP Reagent, Dura**

## 13 . Disposal considerations

- Waste disposal** : The generation of waste should be avoided or minimized wherever possible. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.
- Hazardous waste** : Within the present knowledge of the supplier, this product is not regarded as hazardous waste, as defined by EU Directive 91/689/EEC.

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.

## 14 . Transport information

Regulatory information	UN number	Proper shipping name	Classes	PG*
DOT Classification	Not regulated.	-	-	-
IATA-DGR Class	Not regulated.	-	-	-

PG\* : Packing group

## 15 . Regulatory information

### United States

#### HCS Classification

- : Not regulated.

#### U.S. Federal regulations

- : TSCA 8(a) IUR: sucrose  
**United States inventory (TSCA 8b)**: Not determined.  
 TSCA 12(b) one-time export: 3(2h)-isothiazolone, 2-methyl-  
**SARA 302/304/311/312 extremely hazardous substances**: No products were found.  
**SARA 302/304 emergency planning and notification**: No products were found.  
**SARA 302/304/311/312 hazardous chemicals**: No products were found.  
**SARA 311/312 MSDS distribution - chemical inventory - hazard identification**: No products were found.  
**Clean Water Act (CWA) 307**: No products were found.  
**Clean Water Act (CWA) 311**: No products were found.  
**Clean Air Act (CAA) 112 accidental release prevention**: No products were found.  
**Clean Air Act (CAA) 112 regulated flammable substances**: No products were found.  
**Clean Air Act (CAA) 112 regulated toxic substances**: No products were found.

### Canada

#### WHMIS (Canada)

- : Not controlled under WHMIS (Canada).

#### Canadian lists

- : **CEPA Toxic substances**: None of the components are listed.  
**Canadian ARET**: None of the components are listed.  
**Canadian NPRI**: None of the components are listed.  
**Alberta Designated Substances**: None of the components are listed.  
**Ontario Designated Substances**: None of the components are listed.  
**Quebec Designated Substances**: None of the components are listed.

**Mouse Optimized HRP Reagent, Dura**

**15 . Regulatory information**

**Canada inventory** : Canada inventory: Not determined.

**EU regulations**

**Risk phrases** : This product is not classified according to EU legislation.  
:

**International regulations**

**International lists** : **Australia inventory (AICS)**: All components are listed or exempted.  
**China inventory (IECSC)**: All components are listed or exempted.  
**Korea inventory (KECI)**: Not determined.  
**Philippines inventory (PICCS)**: Not determined.  
**Japan inventory (ENCS)**: Not determined.

**16 . Other information**

**Label requirements** : NOT EXPECTED TO PRODUCE SIGNIFICANT ADVERSE HEALTH EFFECTS WHEN THE RECOMMENDED INSTRUCTIONS FOR USE ARE FOLLOWED.

**Hazardous Material Information System (U.S.A.)** :

Health	0
Flammability	0
Physical hazards	0

The customer is responsible for determining the PPE code for this material.

**National Fire Protection Association (U.S.A.)** :



**Date of printing** : 10/29/2009.

**Date of issue** : 10/29/2009.

**Date of previous issue** : No previous validation.

**Version** : 1

Indicates information that has changed from previously issued version.

**Notice to reader**

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.