



**Item Number** S1498  
**Item** Sodium Thiosulfate, Anhydrous, Reagent  
**CAS Number** 7772-98-7  
**Molecular Formula**  $\text{Na}_2\text{S}_2\text{O}_3$   
**Molecular Weight** 158.11  
Assay..... Min. 97.0%

MAXIMUM LIMITS

Calcium (Ca)..... 0.01%  
Heavy Metals (as Pb)..... 0.001%  
Insoluble Matter..... 0.4%  
Sulfide (S)..... 0.0002%

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:** 14422 S. San Pedro St.  
Gardena, CA 90248  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843

**East Coast Facility:** 755 Jersey Ave.  
New Brunswick, NJ 08901