

# Certificate of Quality



Part Number	<b>MFLX40614-61</b>
Lot Number	<b>SAMPLE</b>
Part Description	<b>FITTING B 1/4X3/16 PVDF</b>
Expiration Date	
Country of Origin	<b>USA</b>

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: Kynar 740
- Material meets USP Class VI, and 21 CFR 177.2510 requirements.
- No animal derived ingredients are used in the manufacture of this product.
- This product does not contain SVHC Candidate list Annex XIV materials as updated July 8th, 2021 (219 substances).
- This product is RoHs compliant - 2002/95/EC without exemptions, 2011/65/EU without exemptions and the latest directive 2015/863 without exemptions.
- Shelf life: 10 years

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.