

# Certificate of Quality



Part Number	MFLX06605-27
Lot Number	0423783232
Part Description	TUBING PTFE 1/16X1/8" 25'
Country of Origin	USA

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: Polytetra-fluoroethylene
- All ingredients used in these products are approved under the FDA CFR 21, subchapter B Part 177.1550 - for use with food and food products.
- The tubing itself has not been tested for approval by the FDA or any other testing body.
- Material is manufactured in an AS9100D with ISO 9001-2015 approved facility.
- All PTFE tubing has been manufactured with virgin resin.
- The proper cleaning of this product for food applications is the responsibility of the end-user.
- Material can be sterilized in an autoclave.
- Material is RoHs Compliant.
- This product is made with a Class VI resin.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.