

# Eppendorf Certificate

## Certificate of Purity

**Product:** Eppendorf Serological Pipet, 10 mL

**Order-No.:** 0030 127.722

**Lot:** L000234G

Products are manufactured in an class 8 clean room which meets the ISO 14644-1 specifications. Testing of the above described lot showed conformity within the limits of detection:

- |             |              |  |
|-------------|--------------|--|
| • Human DNA | Method:      | <i>real-time</i> PCR   |
|             | Limit:       | < 2 pg; less than one human cell   |
| • DNase     | Method:      | DNA Digestion  |
|             | Limit:       | not detectable<br>(LOD 1.0 x 10 <sup>-6</sup> Kunitz units)                                    |
| • RNase     | Method:      | RNA Digestion  |
|             | Limit:       | not detectable<br>(LOD 1.0 x 10 <sup>-6</sup> Kunitz units)                                    |
| • Endotoxin | Method:      | Kinetic-turbidimetric LAL test (Ph. Eur.2.6.14)  |
|             | Limit:       | < 0.01 EU/mL   |
| • Sterility | Method:      | Sterilization with SAL 10 <sup>-6</sup> by irradiation<br>according DIN EN ISO 11137-1:2015-11 |
|             | Test method: | USP, Ph. Eur. 8.0/ 2.6.12  |

Hamburg, January 2022



Joana Tziolis  
Product Life Cycle Manager  
Division Consumables



Monika Schneider  
Vice President Global Quality Management &  
Regulatory Affairs

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

ISO 9001  
Certified

ISO 13485  
Certified

ISO 14001  
Certified