

Rabbit Polyclonal Antibody to P504S α -Methylacyl-CoA Racemase (AMACR)

Catalog No.:	RP 134, RP 134-05
Applications:	Analyte Specific Reagent. Analytical and Performance characteristics are not established. Optimal working conditions must be determined by the end user.
Immunogen:	A synthetic human AMACR peptide.
Host:	Rabbit
Format:	Purified immunoglobulin fraction of rabbit antiserum against AMACR containing sodium azide as a preservative.
Specificity:	P504S has been recently described as a prostate cancer-specific gene that encodes a protein involved in the beta-oxidation of branched chain fatty acids. It stains premalignant lesions of prostate: high-grade prostatic intraepithelial neoplasia (PIN) and atypical adenomatous hyperplasia.
Positive Control:	Prostate carcinoma
Cellular Localization:	Cytoplasmic
Storage:	Store at 2-8°C. Do not use beyond the expiration date stated on the label.
Reference:	Beach et al. Am J Surg Pathol 26 (12): 1588, 2002

ASR: Analyte Specific Reagent

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com