

Product Specification

Test	Specification
Assay (C ₆ H ₁₀ O) (by GC)	≥ 99.0 %
Color (APHA)	≤ 10
Density (g/mL) at 25°C	0.940 – 0.943
Residue after Evaporation	≤ 0.0500 %
Acidity (meq/g)	≤ 0.003
Water (H ₂ O)(by Karl Fischer titrn)	≤ 0.1 %

For Laboratory, Research, or Manufacturing Use

NOTE: Color reported is at time of packaging; material will darken with age.

Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality