



Material	BDH7999-10
Material description	Isopropyl alcohol 70% v/v
Grade	Laboratory Reagent
CAS Number	67-63-0
Molecular formula	CH ₃ CHOHCH ₃
Molecular mass	60.10
Additional information	

Characteristics	Specifications
Assay (V/V)	69.0 --> 71.0 %
Appearance	Clear liquid
Specific gravity	0.873 --> 0.883
Prepared to formulation on file	Confirmed

Signature	Additional information
<p>We certify that this batch conforms to the specifications listed above.</p> <p>This document has been produced electronically and is valid without a signature.</p> <p>Michel Boisvert, Head of the Laboratory - Lachine VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3, QC, Canada</p>	<p>The weights and/or volumes used to prepare this product are NIST traceable. All balances used in the preparation of the product are calibrated weekly against NIST traceable weights. The balances are maintained and serviced on a regular basis by an outside certified company. All volumetric glassware used is NIST traceable and certified as meeting Class A specifications. Unless otherwise indicated above, all chemicals used in the preparation of this product are Reagent ACS grade.</p>

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.