


## Anti-CD45 Antibody [ABT-CD45] (A100028)


### Specifications:

Name:	Anti-CD45 Antibody [ABT-CD45]
Description:	Mouse monoclonal (ABT-CD45) antibody to CD45.
Applications:	IHC
Reactivity:	Human
Immunogen:	Synthetic peptide derived from human CD45.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	ABT-CD45
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Immunogen affinity purified.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline with 0.5% Bovine Serum Albumin and 0.1% Kathon CG/ICP.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

### Images:



Immunohistochemical analysis of paraffin-embedded human Liver using Anti-CD45 Antibody.



Immunohistochemical analysis of paraffin-embedded human Tonsil using Anti-CD45 Antibody.