



Item Number TR106
Item Triacetin, USP
CAS Number 102-76-1
Molecular Formula C₉H₁₄O₆
Molecular Weight 218.21

Assay (Anhydrous Basis).....97.0-100.5%
Specific Gravity @ 25°C.....1.152-1.158
Refractive Index @ 25°C..... 1.429-1.430

MAXIMUM LIMITS

Acidity..... To pass test
Water..... 0.2%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843