

Exon	Mutation	Nucleotide Change	Amino Acid Change	COSMIC ID
Exon 9	E542K	1624G>A	Glutamate (E) to Lysine (K)	760
	E542V	1625A>T	Glutamate (E) to Valine (V)	762
	E542Q	1624G>C	Glutamate (E) to Glutamine (Q)	17442
	E542A	1625A>C	Glutamate (E) to Alanine (A)	1041495
	E542G	1625A>G	Glutamate (E) to Glycine (G)	761
	E545K	1633G>A	Glutamate (E) to Lysine (K)	763
	E545A	1634A>C	Glutamate (E) to Alanine (A)	12458
	E545G	1634A>G	Glutamate (E) to Glycine (G)	764
	E545Q	1633G>C	Glutamate (E) to Glutamine (Q)	27133
	E545D	1635G>T	Glutamate (E) to Aspartate (D)	765
Exon 20	H1047L	3140A>T	Histidine (H) to Lysine (L)	776
	H1047R	3140A>G	Histidine (H) to Arginine (R)	775
	H1047Y	3140A>G	Histidine (H) to Tyrosine (Y)	774