

## VWR® BOD REFRIGERATED INCUBATORS

*Reliable and powerful  
compressor technology*

*Fast cool-down and  
recovery times*

*Broad temperature range  
for application flexibility*





## Excellent stability and value for a wide range of low temperature or BOD applications

### Units Feature:

- Powerful and reliable compressor technology
- Fast cool-down and recovery times, even in hot environments with large sample volume
- Broad temperature range for application flexibility
- Microprocessor control
- Temperature setting sensitivity:  $\pm 0.1^{\circ}\text{C}$
- Push-button controls for temperature setpoint selection
- Forced-air circulation technology
- RTD temperature probe and protected setpoint mode prevent accidental temperature change
- Door key lock protects samples from unauthorized access
- Easy-to-clean, corrosion-resistant construction
- cCSAus listed



### 173 L BOD Refrigerated Incubator

- Economical undercounter design
- High and low safety thermostat backups
- Four adjustable leveling feet for stable setup
- 2A outlet easily supports apparatus inside the unit

### 566 L BOD Refrigerated Incubator

- Ideal for BOD applications and temperature settings at or below ambient
- Temperature setpoint selection with high- and low-temperature protection and simple calibration
- Safety relay and alarm LED alert to over/under temperature conditions
- Cooling switch saves energy at temperatures above ambient\*



### Cat. Number

### Description

89508-420	Refrigerated Incubator 173 L / 6.1 cu ft, 115V 60HZ, temperature range $-10^{\circ}\text{C}$ to $+60^{\circ}\text{C}$
10753-894	Refrigerated Incubator 566 L / 20 cu ft, 115V 60HZ, temperature range $-10^{\circ}\text{C}$ to $+50^{\circ}\text{C}$

\* Based on internal testing

To order, contact your VWR Sales Representative, visit [VWR.com](http://VWR.com), or call 1.800.932.5000 today.



1.800.932.5000 | [vwr.com](http://vwr.com)



Prices and product details are current when published, subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, valid where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered TM or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owners. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.