



FEEL THE SW  
DIFFERENCE.®

# COMFORT AND DEXTERITY FOR EVERYDAY HAND PROTECTION.




TF-95LG

POWDER-FREE NITRILE EXAMINATION GLOVES  
WITH ENERGEL REJUVINATING TECHNOLOGY

**FEATURING AN ENERGEL INNER COATING,  
THESE GLOVES SOOTHE AND MOISTURIZE DRY  
HANDS WITHOUT SACRIFICING DEXTERITY**

- EnerGel® Hand Conditioning Technology clinically proven to hydrate and soothe skin
- Exceeds 1.5 AQL laboratory standards
- Fully textured surface for excellent wet and dry grip
- Touch-screen compatible
- Tear-resistant beaded cuff eases donning
- SHA Dermatologically Approved



Optimal  
thickness  
for durability  
and tactile  
sensitivity

Fully textured  
surface for wet  
and dry grip

Tear-resistant  
beaded cuff for  
easy donning

**TF** **TF-95LG**

WEIGHT	LENGTH	THICKNESS	FEATURES	COLOR
5.0g (large size)	9.5"	3.6 mil	EnerGel®	LIGHT GREEN



# TECHNICAL DATA



## TF-95LG POWDER-FREE NITRILE EXAMINATION GLOVES WITH ENERGEL REJUVINATING TECHNOLOGY

### PRODUCT DESCRIPTION

#### SPECIFICATIONS

PRODUCT	TrueForm®
GRADE	Powder-Free Exam
MATERIAL	Nitrile
THICKNESS	FINGER 4.0mil 0.10mm
	PALM 3.6mil 0.09mm
	CUFF 3.0mil 0.08mm
COLOR	Light Green
SURFACE	OUTER Fully Textured
	INNER EnerGel®
LENGTH	9.5" 242mm
WEIGHT	5.0g
FIT TYPE	Ambidextrous
CUFF TYPE	Beaded
WATERTIGHT AQL	1.5
ORIGIN COUNTRY	China

#### PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	22MPA	20 MPA
ELONGATION:	560%	550%

#### INTEGRATED TECHNOLOGIES

**ENERGEL** Hand Conditioning Technology hydrates and soothes the skin with purified aloe and has been clinically proven to prevent many skin ailments, like dermatitis. It has a healing effect that improves skin health.

Clean Washed  
100% Air Inspected



\* Specifications based on size Large glove

### STANDARDS & CERTIFICATIONS

EN ISO 21420:2020  
EN 455-2000-2015  
ASTM D6319-19  
ASTM F1671-13  
ASTM F739-12  
Reach EU 1907/2006  
EU 2016/425

CAT III Complex Design  
UL NFPA 1999-2018  
FDA 21 CFR 177.2600



Please visit [swsafety.com/symbols](http://swsafety.com/symbols) for information on standards used on this document.

### ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
N128401	XS 5½ – 6	3.0" 75mm	100 gloves per box 10 boxes per case
N128402	S 6½ – 7	3.4" 86mm	100 gloves per box 10 boxes per case
N128403	M 7½ – 8	3.7" 95mm	100 gloves per box 10 boxes per case
N128404	L 8½ – 9	4.2" 107mm	100 gloves per box 10 boxes per case
N128405	XL 9½ – 10	4.6" 117mm	100 gloves per box 10 boxes per case

### APPLICATIONS

#### INDUSTRIAL

Assembly  
Liquid Handling

#### MEDICAL

Chemical Handling  
Clinical  
First Responder EMS  
Handling and Packaging Class 1  
Medical Devices  
Laboratory  
Light Chemical Mixing  
Pharma Production

Keep out of sunlight. Store in a cool, dry place.  
Keep away from sources of ozone or ignition.



R\_N12840X1121