



Item Number P1452  
Item Propylparaben, FCC  
CAS Number 94-13-3  
Molecular Formula  $C_{10}H_{12}O_3$   
Molecular Weight 180.20

Assay ( $C_{10}H_{12}O_3$ ; Dried Basis)..... 99.0-100.5%  
Acidity..... To pass test  
Melting Range..... 95°-98°C

MAXIMUM LIMITS

Lead..... 2 mg/kg  
Loss on Drying..... 0.5%  
Residue on Ignition..... 0.05%

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843