

# ARG49ML



83.75" x 55.25" x 31" (H x W x D)

Performance series Pharma-Lab 49 cu.ft. all-refrigerator in stainless steel with glass doors

### Highlights:

Intelligent microprocessor digital temperature controller

Improved temperature accuracy with a separate temperature probe for continuous product readout

Additional sensor reads interior air temperature for superior temperature control

### Product Features:

|   |   |
|---|---|
| Purpose-built design                            | Designed and engineered to store vaccines and other sensitive materials under stable temperature conditions   |
| Glass doors                                     | Full display of the interior, with an anti-condensate door frame heater to ensure a clear view  |
| Adjustable temperature                          | Interior is pre-set to 5°C, with an adjustable control range from 2 to 10°C   |
| Temperature alarm                               | Visual and audible alarms alert if the interior temperature rises or falls out of set range   |
| Complete stainless steel construction           | Interior and exterior are entirely constructed from stainless steel for lasting durability and easy cleanup   |
| Forced air cooling                              | Top-mounted commercial fan ensures rapid temperature recovery   |
| Automatic defrost with adjustable cycle defrost | Adjustable cycle defrost operation ensures the unit avoid temperature spikes and ice buildup  |
| Self-closing doors                              | Two doors close to save energy but will remain open when the door is opened at a 90° angle for easier loading   |
| Open door alarm                                 | Alarm will beep 3 times if the door is left open for 30 seconds, 5 times at 60 seconds, and will sound continuously if the door is open for more than 5 minutes |

## ARG49ML Specifications:

| Overview                   |                      |
|----------------------------|----------------------|
| Height of Cabinet          | 83.75" (213 cm)      |
| Width                      | 55.25" (140 cm)      |
| Depth                      | 31.0" (79 cm)        |
| Depth with Handle          | 32.75" (83 cm)       |
| Depth with door at 90°     | 56.5" (144 cm)       |
| Capacity                   | 49.0 cu.ft. (1388 L) |
| Defrost Type               | Automatic            |
| Door                       | Glass                |
| Cabinet                    | Stainless Steel      |
| US Electrical Safety       | UL                   |
| Canadian Electrical Safety | ULC                  |
| Sanitation                 | ETL-S                |
| Amps                       | 7.5                  |
| Voltage/Frequency          | 115 V AC/60 Hz       |
| Weight                     | 410.0 lbs. (186 kg)  |
| Shipping Weight            | 430.0 lbs. (195 kg)  |
| Estimated Time to Ship     | 1 to 2 Days          |
| Parts & Labor Warranty     | 2 Years              |
| Compressor Warranty        | 5 Years              |
| Refrigerator Features      |                      |
| Door Swing                 | French               |
| Shelf Type                 | Wire                 |
| Shelf Qty                  | 8                    |
| Adjustable Shelves         | Yes                  |
| Thermostat Type            | Digital              |
| Fan Type                   | Interior             |
| Refrigerant Type           | R290a                |
| Refrigerant Amount         | 3.2oz.               |
| High Side PSI              | 290.0                |
| Low Side PSI               | 155.0                |
| Wheel Qty                  | 4                    |
| Interior Light             | Yes                  |
| Temperature Range          | 2 to 10°C            |
| Dimensions                 |                      |
| Interior Height            | 51.5" (131 cm)       |
| Interior Width             | 50.5" (128 cm)       |
| Interior Depth             | 27.0" (69 cm)        |

