

# Certificate of Quality



Part Number	<b>MFLX31515-89</b>
Lot Number	<b>338313</b>
Part Description	<b>SAN ADAPT.375ID.75 ADCF 10PK</b>
Expiration Date	<b>01FEB2032</b>
Country of Origin	<b>USA</b>

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: P5M6K-080
- The copolymer in this product is FDA 21 CFR 177.1520 (c) 3.1a & has been tested in accordance to USP Class VI and 10993-5 standards.
- This product is not recommended for autoclave sterilization processes in the presence of aqueous solutions.
- Shelf life of this product is recommended at 10 years from the date of manufacturer, if stored in a climate controlled environment out of direct sunlight.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.