

## TECHNICAL DATA SHEET

### Digital Thermometer – Soft Tip

#### Intended Use

RightCare™ Digital Thermometer - Soft Tip, can be used for, safe and accurate measurement of body temperature measurement. Easy to use, and the measurement result is presented on a clear LCD display. The thermometer may be used orally, rectally or under the armpit.

Digital Thermometer – Soft Tip		
Characteristics	Results	
Display resolution	0.1°C	
Measurement Range	89.6~109.4°F	
Accuracy	±0.2°F	
Display Mode	LCD	
Measuring time	1 minute	
Working temperature	10~35°C	
Battery	LR41 type DC 1.5V	
Working Humidity	<80%	
Characteristics	Specification	
Description	Digital Thermometer with soft tip for easy and accurate temperature measurement. Its flexible tip ensures safe temperature measurement. Easy to use, and the measurement result is presented on a clear LCD display.	
Dimensions	Size	N/A

**TECHNICAL DATA SHEET**  
**Digital Thermometer – Soft Tip**



**Visuals**

## TECHNICAL DATA SHEET

### Digital Thermometer – Soft Tip

#### Additional Information

Storage Condition	Factory Certifications	Country of Origin
Store under cool, clean, and dry conditions. Avoid excessive heat (over 55°C)	FDA / CE	China

Packaging	Qty	Dimension	Weight
Primary Packing	320 pcs / case	N/A	6.7 kg
Secondary Packing		15x4x1.5cm	

#### List of Item Codes and Attributes

Item Code / Part Number	Description	Packaging Information	Color	Specific Attributes	Treatment
RCDTST	Digital Thermometer - Soft Tip	320 pcs / case.			

#### Disclaimer :

This document and the information contained herein is the property of RightCare™ and may not be used or reproduced without permission. Product users should conduct all appropriate testing or other evaluations to determine the suitability of RightCare products for a particular purpose or use within a particular environment. RIGHTCARE DISCLAIMS ALL WARRANTIES OTHER THAN WHAT HAS BEEN EXPRESSLY PROVIDED.