

Apr. 9th, 2019

Certificate of Conformity

Product: 0.5 ML CRYOBANK VIAL
Catalogue No.: 374080
Lot No.: 60104948
Date of manufacture: March 2019
Date of expiry: March 2024

Thermo Fisher Scientific confirms that the product supplied has been inspected and tested in accordance to relevant specifications

Test Description	Requirements	Result
Visual	According to specification	Complies
Test of function	According to specification	Complies
USP <85> Bacterial endotoxins test	<0.25 endotoxins unit/ml	Complies
USP <87> "Biological Reactivity test, In Vitro	USP current edition	Complies**
USP<88> "Biological Reactivity test In Vivo, USP Class VI Test"	USP current edition	Complies**
Leak Test	IATA current edition	Complies***
DNase/RNase test	Negative	Complies

**Only performed once

***Performed at certain frequency

Sterility is obtained through irradiation according to ISO 11137 "sterilization of health care products –Radiation" with SAL 10⁻⁶.

The product is CE marked according to EU directive, 2007/47/EC amending 93/42/EC and 98/79/EC for IVD, Medical device.

Quality assurance of the product is performed in accordance with the requirements of the Quality Management System ISO 13485:2016, "Medical devices – Quality management systems – System requirements for regulatory purposes".

Suzhou Site



Jing Zhou
Regulatory Affair Engineer