

Protide Pharmaceuticals, Inc. is devoted to the discovery and development of technologies and products in cell therapy, regenerative medicine, stem cell biology, and cell engineering. Protide has over 35 years of experience in R&D and manufacturing.

TM-235™

Is a serum replacement with similar properties to TCM™ but contains additional components necessary for cells requiring higher concentrations of serum. cGMP Manufactured. FOR CELL CULTURE USE. FOR FURTHER MANUFACTURING USE. NOT FOR HUMAN USE.



Product Highlights	
• Concentration: 50X	• Total Protein: 770ug/ml at 1X Concentrate
• Endotoxin: <0.5EU/ml at 1X	• Storage: 2-8°C
• Protein Origin: Bovine*	• Stability: Minimum 6 months at 2-8°C

*Bovine serum albumin and transferrin are highly purified and heat treated. These components are sourced from documented cattle born, raised and bled in the USA and inspected by the USDA.

VWR Cat. No.	Supplier Cat. No.	Product Volume	Final Packaging Configuration
TM-235™			
76443-258	2001	40ml	Media Bottle
76443-260	2002	100ml	Media Bottle
76443-262	2005	500ml	Media Bottle