

Material Safety Data Sheet

C&S Medium

ACC# 89655

Section 1 - Chemical Product and Company Identification

MSDS Name: C&S Medium**Catalog Numbers:** 2300543**Synonyms:** None.**Company Identification:**

Fisher Diagnostics
 Fisher Scientific Company, LLC
 8365 Valley Pike
 Middletown, VA 22645-0307

For information, call: 800-524-0294**Emergency Number:** 800-524-0294**For CHEMTREC assistance, call:** 800-424-9300**For International CHEMTREC assistance, call:** 703-527-3887

Section 2 - Composition, Information on Ingredients

CAS#	Chemical Name	Percent	EINECS/ELINCS
0-00-0	Agar	0.15	unlisted
367-51-1	Sodium thioglycollate	0.15	206-696-4
7558-79-4	Sodium phosphate dibasic	0.11	231-448-7
143-74-8	Phenol red	0.0024	205-609-7
7732-18-5	Water	Balance	231-791-2

Section 3 - Hazards Identification

EMERGENCY OVERVIEW

Appearance: red gelatinous solution. **Caution!** May cause irritation.**Target Organs:** None.**Potential Health Effects****Eye:** May cause eye irritation.**Skin:** May cause skin irritation.**Ingestion:** May cause irritation of the digestive tract.**Inhalation:** May cause respiratory tract irritation.**Chronic:** No information found.

Section 4 - First Aid Measures