



Item Number BR154
Item Bromothymol Blue, Reagent, ACS
CAS Number 76-59-5
Molecular Formula $C_{27}H_{28}O_5SBr_2$
Molecular Weight 624.39
Clarity of Solution..... To pass test
Visual Transition Interval..... pH 6.0 (yellow) to pH 7.6 (blue)

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843