

Product Name:	PF041WWW/0M
Description:	<p>Corepoint Scientific™ General Purpose Undercounter Freezer Freestanding 4 Cu. Ft. Capacity  One solid swing door, right hinged, not reversible  Mechanical thermostat  Adjustable operating temperature range: -15°C – -25°C  2 fixed shelves  Keyed door lock  Magnetic door gasket for positive seal  Recessed door handle  Leveling legs  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging One year parts and labor warranty  Exterior dimensions: 20 1/4"W x 21 5/8"D x 31 1/2"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost  115V, 60 Hz, .90 Amps, 1/7 HP  UL/C-UL listed  Shipping weight: 94 lbs.</p>