

INCUBATOR-GENIE™

It's a Hybridization Oven, an Incubator, a Rocker/Shaker, a Rotator, ALL-IN-ONE!

75.0°C
28.0°C



Accessories Supplied with Unit



U.S. Pat. 5,918,979

- Ideal for hybridization, digestions, or incubation
- Magnetic stainless steel platform may be easily set to rock or rotate almost any container
- Containers can be easily mounted on the magnetic platform with magnetic clips or strips
- Chamber air is gently and continuously re-circulated ensuring temperature consistency and uniformity

Supplied

Specifications

1 Clip Plate each for 10-13 mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray, Shelf

Temperature

28°C - 75°C (82°F - 167°F)
Accuracy: ± 0.2°C, Uniformity: ± 0.5°C
Resolution: 1.0°C, Stability: ± 0.1°C

Rotating

3-35 RPM

Rocking

6-70 cycles per minute

Platform Capacity

4.5 Kg (10 lb)

Chamber Volume

24 L (0.86 ft³)

Weight

19.5 Kg (43 lb)

Dimensions (D x W x H)

400 x 560 x 370mm

(Overall)

15.75 x 22 x 14.5in

Dimensions (D x W x H)

254 x 362 x 260mm

(Chamber)

10 x 14.25 x 10.25in

Alarms

Visible and Audible for Temperature

Warranty

2 Years



Fig.	VWR Cat. No.	Description	Fig.	VWR Cat. No.	Description
--	97060-138	Incubator Genie, 120V	8	15001-016	Bag Expansion Kit
1	15001-008	Clip Plate for 12 each of 10-13mm tubes, 1/pk	9	14003-576	Clip Plates for Hybridization Tubes, 2/pk
2	15001-010	Clip Plate for 6 each of 15-17mm tubes, 1/pk	10	14003-582	Stackable Shelf, 1/pk
3	15001-012	Clip Plate for 3 each of 28-30mm tubes, 1/pk	11	14222-140	Universal Clip Plates, 2/pk
4	15001-002	Clip Plates for 12 each of 10-13mm tubes, 6/pk	12	14222-144	Magnetic Attachable Covered 100 Microtube Box, 1/pk
5	15001-004	Clip Plates for 6 each of 15-17mm tubes, 6/pk	13	14003-574	3-D Orbital Shaker Attachment
6	15001-006	Clip Plates for 3 each of 28-30mm tubes, 6/pk	14	490015-482	Horizontal Orbital Shaker Attachment
7	15001-014	Bag Mounting Strips, 16/pk			