



Item Number G1040
Item Glycerin, Synthetic, FCC
CAS Number 56-81-5
Molecular Formula $C_3H_8O_3$
Molecular Weight 92.10

Assay ($C_3H_8O_3$)..... 95.0-100.5%
Specific Gravity @ 25°C..... Min. 1.249
Color..... To pass test

MAXIMUM LIMITS

Chlorinated Compounds (as Cl)..... 0.003%
Fatty Acids and Esters (as Butyric Acid; Limit about 0.1%)..... To pass test
Lead..... 1 mg/kg
Readily Carbonizable Substances..... To pass test
Residue on Ignition..... 0.01%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843