

Product Name:	NSWDR492SSS/0
Description:	<p>Corepoint Scientific™ White Diamond Series Laboratory and Medical Stainless Steel Refrigerator 49 Cu. Ft. capacity Two swing solid doors, self-closing 10 ea. Adjustable Shelves Microprocessor temperature controller Touch Screen Digital Temperature Display with user friendly interface Alarm Module with battery back-up to support temperature monitoring and alarms Password-protected control panel C°/F° switchable 0.1° resolution Adjustable operating temperature: 1°C to 10°C (34°F to 50°F) Active temperature scale bar display Resettable 1-point min/max history Digital 24-hour temperature chart Temperature data logging USB port for data transfer and software updates Real-time operational display clock 4-20mA output Adjustable alarm volume Audible and visual high/low temperature alarms 2 probes; 1 in air and 1 in sample bottle Continuous product (simulated) temperature display Sensor error alarm Power failure alarm Door ajar alarm Low battery alarm Alarm validation (audible, visual, and remote alarm) Alarm mute with ring back Remote alarm contacts Customizable alarm time delay settings Forced draft air circulation Keyed door locks Magnetic door gaskets for positive seal Integrated door handles LED interior lights are safety shielded and switch controlled Probe access port (3/4") Stainless steel exterior White-powder coated aluminum interior with 304 stainless steel floor Easily cleaned smooth interior Hydrocarbon natural refrigerant (R290) Three year parts and labor warranty, plus an additional two year compressor parts warranty 4 swivel casters (5 3/4") Exterior dimensions: 54" W x 33 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ETL listed Shipping weight: 479 lbs.</p>