



Hexanes AR (ACS)

Product No. 5189
Specifications current as of: Jun 8 2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Assay (sum of 5 isomers, total hexanes, plus methylcyclopentane)(by GC)	min. 98.5 %
Assay (as n-hexane) (by GC, corrected for water)	Actual Value Reported
Appearance (clear, colorless liquid)	Passes Test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001 %
Sulfur Compounds (as S)	max. 0.005 %
Thiophene	Passes Test
Water-Soluble Titrable Acid (meq/g)	max. 0.0003

For Laboratory, Research or Manufacturing Use