

Catalog Number: BDH20067.400

Catalog Description: Acetone

Grade: HiPerSolv CHROMANORM for HPLC

Date Current Specification Set: 2-16-16


All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics

Specifications

Assay (by GC-FID) (Corrected for Water)	Min 99.5 %
Acidity	Max 0.0005 meq/g
Alkalinity	Max 0.0002 meq/g
Residue after evaporation Water	Max 5 ppm
(by KF coulometry)	Max 0.5 %
Transmittance at	
345 nm	Min 50.0 %
350 nm	Min 80.0 %
360 nm	Min 98.0 %

Signed on behalf of VWR: _____

Quality Control Manager 

Additional Information (optional):