

Citric Acid 50%, Solution
Biotech Reagent



Material No.: 0341-27

Revision No.: 0

Product Specification

, GMP Manufactured Product

Test	Specification
Assay (HOC(COOH)(CH ₂ COOH) ₂)	48 - 52 %
Identification	Passes Test
ACS - Chloride (Cl)	≤ 0.001 %
Heavy Metals (as Pb)	≤ 0.001 %
ACS - Iron (Fe)	≤ 3 ppm
Phosphate (PO ₄)	≤ 0.001 %

GMP Manufactured Product

CAUTION: For Manufacturing, processing or repackaging

Made from Purified Water U.S.P. and Citric Acid, anhydrous U.S.P./A.C.S.

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Filtered through a 0.2 micron filter.

Packaging Site: Paris Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700