



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Revision Date 12/11/2014

Version 1.1

## SECTION 1. Identification

### Product identifier

Product number	106290
Product name	Sodium benzoate suitable for use as excipient EMPROVE® exp Ph Eur,BP,NF,FCC,E 211
CAS-No.	532-32-1

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Food and pharmaceutical industry, Cosmetic raw material
-----------------	---

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   General Inquiries: +1-978-715-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)
---------	--

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

## SECTION 2. Hazards identification

### GHS-Labeling

Not a dangerous substance according to GHS.

### Other hazards

None known.

## SECTION 3. Composition/information on ingredients

Formula	$C_6H_5COONa$	$C_7H_5NaO_2$ (Hill)
Molar mass	144.10 g/mol	

Remarks	No hazardous ingredients according to the OSHA Hazard Communication Standard 29 CFR 1910.1200.
---------	---

## SECTION 4. First aid measures

### Description of first-aid measures

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

106290

Version 1.1

Product name

Sodium benzoate suitable for use as excipient EMPROVE® exp Ph  
Eur,BP,NF,FCC,E 211

---

### *Inhalation*

After inhalation: fresh air.

### *Skin contact*

In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/  
shower.

### *Eye contact*

After eye contact: rinse out with plenty of water.

### *Ingestion*

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

Never give anything by mouth to an unconscious person.

### **Most important symptoms and effects, both acute and delayed**

Nausea, Vomiting, Stomach/intestinal disorders

### **Indication of any immediate medical attention and special treatment needed**

No information available.

---

## **SECTION 5. Fire-fighting measures**

### **Extinguishing media**

#### *Suitable extinguishing media*

Water, Carbon dioxide (CO<sub>2</sub>), Foam, Dry powder

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### **Special hazards arising from the substance or mixture**

Combustible.

Risk of dust explosion.

Development of hazardous combustion gases or vapors possible in the event of fire.

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

In the event of fire, wear self-contained breathing apparatus.

#### *Further information*

Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Avoid inhalation of dusts. Evacuate the danger area,  
observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### **Environmental precautions**

Do not empty into drains.

### **Methods and materials for containment and cleaning up**

Cover drains. Collect, bind, and pump off spills.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

106290

Version 1.1

Product name

Sodium benzoate suitable for use as excipient EMPROVE® exp Ph  
Eur,BP,NF,FCC,E 211

---

Observe possible material restrictions (see sections 7 and 10).

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

### Conditions for safe storage, including any incompatibilities

Tightly closed. Dry.

Store at +5°C to +30°C (+41°F to +86°F).

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

Contains no substances with occupational exposure limit values.

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### *Hygiene measures*

Change contaminated clothing. Wash hands after working with substance.

### *Eye/face protection*

Safety glasses

### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### *Respiratory protection*

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

## SECTION 9. Physical and chemical properties

Physical state

solid

Color

white

Odor

odorless

Odor Threshold

Not applicable

---

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

106290

Version 1.1

Product name

Sodium benzoate suitable for use as excipient EMPROVE® exp Ph  
Eur,BP,NF,FCC,E 211

---

pH	ca. 9 at 100 g/l 68 °F ( 20 °C)
Melting point	410 - 430 °C
Boiling point/boiling range	Not applicable
Flash point	> 212 °F ( > 100 °C)
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	Not applicable
Upper explosion limit	Not applicable
Vapor pressure	Not applicable
Relative vapor density	Not applicable
Density	1.44 g/cm <sup>3</sup> at 68 °F ( 20 °C)
Relative density	No information available.
Water solubility	660 g/l at 68 °F ( 20 °C)
Partition coefficient: n-octanol/water	log Pow: -2.27 (calculated) (Lit.) Bioaccumulation is not expected.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Oxidizing properties	none
Ignition temperature	> 932 °F ( > 500 °C)
Bulk density	ca. 350 kg/m <sup>3</sup>

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

106290

Version 1.1

Product name

Sodium benzoate suitable for use as excipient EMPROVE® exp Ph  
Eur,BP,NF,FCC,E 211

---

## SECTION 10. Stability and reactivity

### Reactivity

Risk of dust explosion.

### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

### Possibility of hazardous reactions

Violent reactions possible with:

Strong oxidizing agents, Strong acids

### Conditions to avoid

no information available

### Incompatible materials

no information available

### Hazardous decomposition products

no information available

---

## SECTION 11. Toxicological information

### Information on toxicological effects

#### *Likely route of exposure*

Eye contact, Skin contact, Ingestion

#### *Acute oral toxicity*

LD50 Rat: 3,140 mg/kg (External MSDS)

Symptoms: Nausea, Vomiting, Stomach/intestinal disorders

#### *Skin irritation*

Rabbit

Result: No irritation

OECD Test Guideline 404

#### *Eye irritation*

Rabbit

Result: slight irritation

OECD Test Guideline 405

#### *Genotoxicity in vitro*

Ames test

Result: negative

(Lit.)

#### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

#### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

106290

Version 1.1

Product name

Sodium benzoate suitable for use as excipient EMPROVE® exp Ph  
Eur,BP,NF,FCC,E 211

---

## *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

## **Carcinogenicity**

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

NTP

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

ACGIH

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

## **Further information**

Hazardous properties cannot be excluded but are unlikely when the product is handled appropriately.

Handle in accordance with good industrial hygiene and safety practice.

---

## **SECTION 12. Ecological information**

### **Ecotoxicity**

#### *Toxicity to fish*

LC50 *Leuciscus idus* (Golden orfe): 460 mg/l (ECOTOX Database)

#### *Toxicity to daphnia and other aquatic invertebrates*

EC50 *Daphnia*: > 100 mg/l; 48 h (External MSDS)

#### *Toxicity to algae*

IC50 algae: > 10 mg/l; 72 h (External MSDS)

#### *Toxicity to bacteria*

EC0 Bacteria: 1,000 mg/l (External MSDS)

### **Persistence and degradability**

#### *Biodegradability*

90 %; 7 d

OECD Test Guideline 301B

Readily biodegradable.

### **Bioaccumulative potential**

#### *Partition coefficient: n-octanol/water*

log Pow: -2.27

(calculated)

(Lit.) Bioaccumulation is not expected.

### **Mobility in soil**

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

106290

Version 1.1

Product name

Sodium benzoate suitable for use as excipient EMPROVE® exp Ph  
Eur,BP,NF,FCC,E 211

---

No information available.

*Additional ecological information*

Discharge into the environment must be avoided.

---

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

## SECTION 14. Transport information

### Land transport (DOT)

Not classified as dangerous in the meaning of transport regulations.

### Air transport (IATA)

Not classified as dangerous in the meaning of transport regulations.

### Sea transport (IMDG)

Not classified as dangerous in the meaning of transport regulations.

---

## SECTION 15. Regulatory information

### United States of America

#### SARA 313

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

#### SARA 302

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

#### Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

#### DEA List I

Not listed

#### DEA List II

Not listed

### US State Regulations

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

106290

Version 1.1

Product name

Sodium benzoate suitable for use as excipient EMPROVE® exp Ph  
Eur,BP,NF,FCC,E 211

---

## Massachusetts Right To Know

### Remarks

No components are subject to the Massachusetts Right to Know Act.

## California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

## Notification status

TSCA:

All components of the product are listed in the TSCA-inventory.

DSL:

All components of this product are on the Canadian DSL.

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Revision Date 12/11/2014

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*