

MATERIAL SAFETY DATA SHEET

Animal Tissues

Rockland Immunochemicals, Inc.
P.O. Box 326
Gilbertsville, PA 19512

Emergency Telephone Number:
Chemtrec (800) 424-9300

Information Telephone Number:
(800) 656-7625

DATE PREPARED: 08-18-2006

REVIEWED:

SECTION 1:-----CHEMICAL IDENTIFICATION-----

NAME: ANIMAL TISSUES

SECTION 2:-----COMPOSITION/INFORMATION ON INGREDIENTS-----

SUBSTANCE NAME CAS #
ANIMAL TISSUES

SECTION 3:-----HAZARDS IDENTIFICATION-----

HMIS RATING

HEALTH: 0
FLAMMABILITY: 0
REACTIVITY: 0

NFPA RATING

HEALTH: 0
FLAMMABILITY: 0
REACTIVITY: 0

FOR ADDITIONAL INFORMATION ON TOXICITY, PLEASE REFER TO SECTION 11.

SECTION 4:-----FIRST-AID MEASURES-----

ORAL EXPOSURE

IF SWALLOWED, WASH OUT MOUTH WITH WATER PROVIDED PERSON IS CONSCIOUS. CALL A PHYSICIAN IMMEDIATELY.

INHALATION EXPOSURE

IF INHALED, REMOVE TO FRESH AIR. IF NOT BREATHING GIVE ARTIFICIAL RESPIRATION. IF BREATHING IS DIFFICULT, GIVE OXYGEN.

DERMAL EXPOSURE

IN CASE OF SKIN CONTACT, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT LEAST 15 MINUTES. REMOVE CONTAMINATED CLOTHING AND SHOES. CALL A PHYSICIAN.

EYE EXPOSURE

IN CASE OF CONTACT WITH EYES, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT LEAST 15 MINUTES. ASSURE ADEQUATE FLUSHING BY SEPARATING THE EYELIDS WITH FINGERS. CALL A PHYSICIAN.

SECTION 5:-----FIRE FIGHTING MEASURES-----

AUTOIGNITION TEMP

N/A

FLAMMABILITY

N/A

EXTINGUISHING MEDIA

SUITABLE: WATER SPRAY. CARBON DIOXIDE, DRY CHEMICAL POWDER, OR APPROPRIATE FOAM.

FIREFIGHTING

PROTECTIVE EQUIPMENT: WEAR SELF-CONTAINED BREATHING APPARATUS AND PROTECTIVE CLOTHING TO PREVENT CONTACT WITH SKIN AND EYES.
SPECIFIC HAZARD(S): EMITS TOXIC FUMES UNDER FIRE CONDITIONS.

SECTION 6:-----ACCIDENTAL RELEASE MEASURES-----
METHODS FOR CLEANING UP
VENTILATE AREA AND WASH SPILL SITE AFTER MATERIAL PICKUP IS COMPLETE.

SECTION 7:-----HANDLING AND STORAGE-----
HANDLING
USER EXPOSURE: AVOID INHALATION. AVOID CONTACT WITH EYES, SKIN,
AND CLOTHING. AVOID PROLONGED OR REPEATED EXPOSURE.

SECTION 8:-----EXPOSURE CONTROLS/PERSONAL PROTECTION-----
ENGINEERING CONTROLS
SAFETY SHOWER AND EYE BATH.
PERSONAL PROTECTIVE EQUIPMENT
RESPIRATORY: USE RESPIRATORS AND COMPONENTS TESTED AND APPROVED
UNDER APPROPRIATE GOVERNMENT STANDARDS SUCH AS NIOSH (US) OR CEN (EU).
RESPIRATORY PROTECTION IS NOT REQUIRED. WHERE PROTECTION IS
DESIRED, USE MULTI-PURPOSE COMBINATION (US) OR TYPE ABEK (EN 14387) RESPIRATOR
CARTRIDGES.
HAND: PROTECTIVE GLOVES.
EYE: CHEMICAL SAFETY GOGGLES.
GENERAL HYGIENE MEASURES
WASH THOROUGHLY AFTER HANDLING.

SECTION 9:-----PHYSICAL AND CHEMICAL PROPERTIES-----
APPEARANCE AND ODOR
ANIMAL TISSUE

SECTION 10:-----STABILITY AND REACTIVITY -----
STABILITY
STABLE

SECTION 11:-----TOXICOLOGICAL INFORMATION-----
ROUTE OF EXPOSURE
SKIN CONTACT: MAY CAUSE SKIN IRRITATION.
SKIN ABSORPTION: MAY BE HARMFUL IF ABSORBED THROUGH THE SKIN.
EYE CONTACT: MAY CAUSE EYE IRRITATION.
INHALATION: MATERIAL MAY BE IRRITATING TO MUCOUS MEMBRANES AND UPPER
RESPIRATORY TRACT. MAY BE HARMFUL IF INHALED.
INGESTION: MAY BE HARMFUL IF SWALLOWED.
SIGNS AND SYMPTOMS OF EXPOSURE:
TO THE BEST OF OUR KNOWLEDGE, THE CHEMICAL, PHYSICAL, AND TOXICOLOGICAL
PROPERTIES HAVE NOT BEEN THOROUGHLY INVESTIGATED.

SECTION 12:-----ECOLOGICAL INFORMATION-----
DATA NOT YET AVAILABLE

SECTION 13:-----DISPOSAL CONSIDERATIONS-----
OBSERVE ALL FEDERAL, STATE, AND LOCAL ENVIRONMENTAL REGULATIONS.

SECTION 14:-----TRANSPORT INFORMATION-----
DATA NOT YET AVAILABLE

SECTION 15:-----REGULATORY INFORMATION-----
NOT AVAILABLE.

SECTION 16:-----OTHER INFORMATION-----
THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT PURPORT TO BE
ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. MANUFACTURER SHALL NOT BE
HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING FROM CONTACT WITH THE
ABOVE PRODUCT.

CHECK OUT OUR:

[PHOSPHO SPECIFIC ANTIBODIES](#)

[MOLECULAR BIOLOGY BUFFERS](#)

[CUSTOM ANTIBODY SERVICES](#)