

Certificate of Quality



Part Number	MFLX06411-58
Lot Number	28509935
Part Description	TUBING SILICONE 5/8X7/8IN 25FT
Expiration Date	31DEC2026
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: Peroxide Cured Silicone
- Per the manufacturers technical data sheet:
- FDA Approved for food contact
- USP Class VI - OK
- Autoclavable - OK
- This product is not approved or intended for, and should not be used for medical, clinical, surgical or other patient oriented applications.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.