

## Product Specification

Meets ACS Reagent Chemical Requirements,

Test	Specification
Assay (C <sub>6</sub> H <sub>5</sub> CH <sub>3</sub> ) (by GC)	≥ 99.5 %
Color (APHA)	≤ 10
Residue after Evaporation	≤ 0.001 %
Substances Darkened by H <sub>2</sub> SO <sub>4</sub>	Passes Test
Sulfur Compounds (as S)	≤ 0.003 %
Water (by KF, coulometric)	≤ 0.03 %
Benzene (C <sub>6</sub> H <sub>6</sub> ), For Information Only	ppm

For Laboratory, Research, or Manufacturing Use

Date Received: \_\_\_ / \_\_\_ / \_\_\_

Date Opened: \_\_\_ / \_\_\_ / \_\_\_

Meets Reagent Specifications for testing USP/NF monographs

Packaging Site: Paris Mfg Ctr & DC