


Anti-Actin Antibody [ABT-ACTN] (A98880)


Specifications:

Name:	Anti-Actin Antibody [ABT-ACTN]
Description:	Mouse monoclonal (ABT-ACTN) antibody to Actin.
Applications:	IHC, WB
Reactivity:	Human
Immunogen:	Synthetic peptide derived from human Actin pan.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	ABT-ACTN
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Immunogen affinity purified.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline with 0.5% Bovine Serum Albumin and 0.1% Kathon CG/ICP.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:




Immunohistochemical analysis of paraffin-embedded human Cardiac muscle using Anti-Actin Antibody.




Immunohistochemical analysis of paraffin-embedded human Colon using Anti-Actin Antibody.

Anti-Actin Antibody [ABT-ACTN] (A98880)

Images continued:



Immunohistochemical analysis of paraffin-embedded human Colon using Anti-Actin Antibody.



Immunohistochemical analysis of paraffin-embedded human Colon using Anti-Actin Antibody.