

Universal-Mount Eco Glove Dispensers

Universal-Mount Eco Glove Dispensers help your employees and the environment.

All Eco-Dispensers are made from .150" 100% recycled PETG providing a more durable and eco-friendly alternative to acrylic dispensers. Single, double and triple dispensers are available in clear and white. Optional magnet mounts allow easy placement for use on filing cabinets, storage cabinets, metal doors, etc. One magnet required for single, two for double and three for triple dispensers. Made in the USA.



Triple-Box Dispenser, White
(76026-304)



Double-Box Dispenser, White
(76026-302)



Single-Box Dispenser, Clear
(76026-178)

- All Dispensers Mount Either Vertically or Horizontally
- Available in Clear & White
- Optional Magnet Mount for Metal Surfaces
Order 1 for single, 2 for double, 3 for triple

VWR Cat. No.	Description	Length	Width	Height
76026-178	Glove Dispenser-Single Box, Clear	10.5"	6.0"	4.0"
76026-180	Glove Dispenser-Double Box, Clear	10.5"	12.0"	4.0"
76026-182	Glove Dispenser-Triple Box, Clear	10.5"	16.0"	4.0"
76026-184	Glove Dispenser-Single Box, White	10.5"	6.0"	4.0"
76026-302	Glove Dispenser-Double Box, White	10.5"	12.0"	4.0"
76026-304	Glove Dispenser-Triple Box, White	10.5"	16.0"	4.0"
76026-306	Optional Magnet Mount	10.0"	5.0"	.125"