



FEEL THE SW  
DIFFERENCE®

# DESIGNED TO ENHANCE WORKER PERFORMANCE.


## PowerForm®

### PF-95TL

POWDER-FREE ECOTEK BIODEGRADABLE  
NITRILE EXAMINATION GLOVES\*

**VERSATILE BIODEGRADABLE 5 MIL SINGLE-USE  
GLOVE COMBINES MAXIMUM COMFORT AND  
CONTROL WITH HIGH-QUALITY CONSTRUCTION**

- EcoTek® Biodegradable Technology accelerates landfill breakdown without any performance loss\*
- Silicone-free for critical environments
- Exceeds 1.5 AQL laboratory standards







Optimal  
thickness for  
durability and  
tactile sensitivity

Fully textured  
surface for wet  
and dry grip

Tear-resistant  
beaded cuff for  
easy donning

**PF** PF-95TL

WEIGHT	LENGTH	THICKNESS	FEATURES	COLOR
 6.5g (large size)	 9.5"	 5.0 mil	EcoTek® Clean Washed Air Inspected	 TEAL

**ECOTEK**



## TECHNICAL DATA



**PF-95TL**  
POWDER-FREE ECOTEK BIODEGRADABLE  
NITRILE EXAMINATION GLOVES\*

### PRODUCT DESCRIPTION

#### SPECIFICATIONS

<b>PRODUCT</b>	PowerForm®	
<b>GRADE</b>	Powder-Free Exam	
<b>MATERIAL</b>	Nitrile	
<b>THICKNESS</b>	<b>FINGER</b>	5.8mil 0.15mm
	<b>PALM</b>	5.0mil 0.13mm
	<b>CUFF</b>	3.9mil 0.10mm
<b>COLOR</b>	Teal	
<b>SURFACE</b>	Fully Textured	
<b>LENGTH</b>	9.5"	242mm
<b>WEIGHT</b>	6.5g	
<b>FIT TYPE</b>	Ambidextrous	
<b>CUFF TYPE</b>	Beaded	
<b>WATERTIGHT AQL</b>	1.5	
<b>ORIGIN COUNTRY</b>	China	

\* Specifications based on size Large glove

#### PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
<b>TENSILE</b>		
<b>STRENGTH:</b>	21 MPA	20 MPA
<b>ELONGATION:</b>	560%	540%

#### INTEGRATED TECHNOLOGIES

**ECOTEK** Biodegradable Technology: 92.6% biodegradable in 2.5 years. ASTM D5526-12 test results do not predict future biodegradation.\*

100% Air Inspected


Clean Washed



### STANDARDS & CERTIFICATIONS

EN ISO 21420:2020  
EN 455-2000-2015  
ASTM D6319-19  
ASTM D6978-05 (2019)  
ASTM F1671-13  
ASTM F739-12  
Reach EU 1907/2006

EU 2016/425  
CAT III Complex Design



Please visit [swsafety.com/symbols](http://swsafety.com/symbols) for information on standards used on this document.

### ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
N200362	S 6½ – 7	3.5" 89mm	100 gloves per box 10 boxes per case
N200363	M 7½ – 8	3.8" 98mm	100 gloves per box 10 boxes per case
N200364	L 8½ – 9	4.3" 108mm	100 gloves per box 10 boxes per case
N200365	XL 9½ – 10	4.6" 116mm	100 gloves per box 10 boxes per case
N200366	2XL 10½ – 11	4.9" 124mm	100 gloves per box 10 boxes per case

### APPLICATIONS

**INDUSTRIAL**  
Food Handling and Packaging  
Janitorial Maintenance  
Refinery/Oil

**MEDICAL**  
Chemical Handling

Keep out of sunlight. Store in a cool, dry place.  
Keep away from sources of ozone or ignition.



R\_N20036X1121



\*Under test method ASTM D5526-12, a standard test method for determining anaerobic biodegradation of plastics in accelerated landfill conditions. Current biodegradation results do not predict future results.