

# The Gibraltar SteriKit®



## BENEFITS

- Offered in a variety of sizes (2-50 cc vials – clear or amber)
- Ready to use
- Final packaging under ISO 5 conditions
- Customizable
- Saves time, space and money
- Shorter lead times
- Sterile vacuum packaging
- Easy to use disposable packaging
- Certified in-house
- No ETO residual
- No trembling during transport
- Made and shipped in the USA

## QUALITY STANDARDS

- Backed by Gibraltar Laboratories, the award winning, 46 year old 2015 ISPE nominated Facility of the Year.
- Every SteriKit® is certified by Gibraltar Laboratories sterile, pyrogen and particulate free (USP <71>, <85> & <788>).
- These quality attributes for the product are important for patient safety.

## Why Gibraltar SteriKit® is 100% The Right Choice:

- Vials, stoppers and seals certified in-house as sterile, pyrogen and particulate free.
- Manufactured and packaged in cGMP ISO clean rooms.
- Reliable delivery time.
- No SteriKit® customer has ever experienced a production delay.
- Full analytical service to support development from pre-clinical R&D through commercial manufacturing and sterilization.
- Processing vials in house is not optimal. Repurpose space to increase production.
- Gibraltar does not compete with Contract Manufacturing Organizations (CMOs).
- No order is too small.

Take control and accelerate your manufacturing cycle with the **Gibraltar Advantage®**:

Sterile vacuum packaging, shorter lead times and more customizable solutions.



For ordering information:

1.800.932.5000  
www.vwr.com

# The Gibraltar SteriKit®

## APPLICATIONS

- Aseptic fill of injectable liquids
- Vaccine companies
- Biopharmaceutical companies
- 503B Compounding pharmacies
- Small batch sizes
- Orphan drugs
- Personalized medicines
- Clinical trials
- Technical batches
- Diagnostics

## EXCEPTIONAL PROCESSING

- Type I borosilicate vials are USP WFI washed, aseptically repackaged in a 316 SS tray and USP dry heat depyrogenated.
- The stoppers and seals are USP WFI washed, aseptically repackaged and cGMP moist heat steam sterilized.
- Following sterilization, the components are assembled, vacuum packed and heat sealed in a pre-sterilized LDPE bag under ISO Class 5 conditions.

## LEAVE IT TO THE PROS

Particulate problems? ETO residual? Endotoxin problems? Lengthy lead times and high outsourcing costs? Autoclave issues? Wasting time, space and money? At Gibraltar Laboratories, Division of Sterilization Services; ***we heard you.***

Don't waste time, space and money washing, depyrogenating, sterilizing and certifying vials, seals and stoppers. We want you to use that time, space and money to do what you do best – accelerate your manufacturing process. Whether it be a diagnostic, injectable drug or an implant, Gibraltar is here to help. Through our world class sterilization services, Gibraltar offers pyrogen free, particulate free and sterile vials, seals and stoppers that are ready to be inserted directly into your manufacturing facility.

Our facilities are modern and our processes are validated and compliant as per cGMP standards. Our packaging is easy to use, clean, disposable and prevents vial trembling during shipping. The clear final packaging also allows for visual inspection of the vials before use.

Give us a call today and ask how we can help with your sterilization needs and how our **SteriKit®** will save you time, space and money™.

## TRUST US

The Gibraltar SteriKit® is a ready-to-use kit that combines sterile, pyrogen and particulate-free vials, stoppers, and seals in sterile vacuum-sealed, disposable packaging. The SteriKit® is highly-customizable and scalable to meet your specific needs ranging from small, one-off R&D projects to large clinical and commercial manufacturing. Since development and manufacturing demands vary by client, product and application; the SteriKit® provides a robust platform that can help reduce variability and minimize risk to your products, patients, and reputation. It just makes good business sense.

**We believe that our client's success is our success.** Take control and accelerate your manufacturing with the Gibraltar SteriKit®

## Ordering Information

Product	Size	VWR Cat. No.
SteriKit® Full Kit	2 mL	75805-150
SteriKit® Full Kit	3 mL	75805-152
SteriKit® Full Kit	5 mL	75805-154
SteriKit® Full Kit	10 mL	75805-156
SteriKit® Full Kit	20 mL	75805-158
SteriKit® Vials & Stoppers Only	2 mL	75805-160
SteriKit® Vials & Stoppers Only	3 mL	75805-178
SteriKit® Vials & Stoppers Only	5 mL	75805-162
SteriKit® Vials & Stoppers Only	10 mL	75805-164
SteriKit® Vials & Stoppers Only	20 mL	75805-166
SteriKit® Vials Only	2 mL	75805-168
SteriKit® Vials Only	3 mL	75805-170
SteriKit® Vials Only	5 mL	75805-172
SteriKit® Vials Only	10 mL	75805-174
SteriKit® Vials Only	20 mL	75805-176

## Full Offering of Steam Sterilization Services

Cycle Development ◦ Validation ◦ Air Over Pressure ◦ Pulse Pre-Vac

...call us today!

