



**Item Number** AD130  
**Item** Adipic Acid, FCC  
**CAS Number** 124-04-9  
**Molecular Formula** C6H10O4  
**Molecular Weight** 146.14

Assay (C<sub>6</sub>H<sub>10</sub>O<sub>4</sub>; Anhydrous Basis)..... 99.6-101.0%  
Melting Range..... 151.5°-154°C

MAXIMUM LIMITS

Lead..... 2 mg/kg  
Residue on Ignition..... 0.002%  
Water..... 0.2%

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843