

## Product Specification

Test	Specification
Assay (Al <sub>2</sub> O <sub>3</sub> )	≥ 98.5 %
Loss on Ignition	≤ 1.0 %
Chloride (Cl)	≤ 0.005 %
Sulfate (SO <sub>4</sub> )	≤ 0.005 %
Heavy Metals (as Pb)	≤ 0.001 %
Iron (Fe)	≤ 0.02 %

For Laboratory, Research, or Manufacturing Use

Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier  
Vice President Global Quality