

A person wearing blue nitrile gloves is using a clear pipette to transfer a liquid into a clear plastic vial. The vial has a white lid that is open. The background is a blurred blue color.

**avantor**<sup>™</sup>  
delivered by **VWR**<sup>™</sup>

[VWR.COM](http://VWR.COM)

**DILU-LOK**<sup>™</sup>

Pre-filled Dilution Vials

# *With Hardy, dilution vials are a snap!*



The Dilu-Lok™ eliminates the time and money spent on preparing your own dilution vials! Our pre-filled dilution vial is designed especially for the food, dairy, pharmaceutical, and water industries. Our easy to use snap lid vials are filled to 25ml, 90ml, or 99ml with a variety of buffering and enrichment solutions. Each lot is tested for pH, sterility, and fill volume.



Convenient Write-on area



Easy-to-remove, shrink wrap

**Heavy duty** - Vials are made of sturdy polypropylene with a leak-resistant seal that holds up even during air transport.

**Large Opening** - A wide 45mm opening accommodates most sample sizes. Each vial is clearly marked with a 90ml and 99ml fill line.

**One-Handed** - Flip-top lid allows for easy one-handed opening and closing to make pipetting easier. Saves unnecessary motion! No awkward screw caps or plastic strips to tear off.

**Variety** - Available with peptone water, TAT broth, buffered peptone water, butterfield's buffer, deionized water, phosphate buffer with MgCl, and as an empty, sterile container.

**Leak-Resistant** - Re-engineered for an improved seal. Perforated shrink wrap is easy to remove and ensures product integrity.

Description	VWR Cat. no.
Butterfield's Phosphate Buffer, 25ml fill, 50/pk	10789-398
Butterfield's Phosphate Buffer, 90ml fill, 50/pk	89407-110
Butterfield's Phosphate Buffer, 99ml fill, 50/pk	89407-112
Peptone Water 0.1%, 90ml fill, irradiated, 50/pk	89407-106
Peptone Water 0.1%, 99ml fill, irradiated, 50/pk	89407-108
Phosphate Buffer with Magnesium Chloride, 90ml fill, 50/pk	89407-114
Phosphate Buffer with Magnesium Chloride, 99ml fill, 50/pk	89407-116
TAT (Tryptone Azolectin Tween®) Broth, 90ml fill, irradiated, 50/pk	10753-958
Buffered Peptone Water, 90ml fill, irradiated, 50/pk	89407-096
Buffered Peptone Water, 99ml fill, irradiated, 50/pk	10145-966
Buffered Peptone Water with 1% Tween® 20, 90ml fill, irradiated, 50/pk	89407-098
Deionized Water, 50ml fill, 50/pk	76169-782
Deionized Sterile Water, 90ml fill, 50/pk	89407-100
Deionized Sterile Water, 99ml fill, 50/pk	89407-102
Dilu-Lok™ container, Empty, Sterile, 250/pk	89423-184







**avantor**<sup>™</sup>  
delivered by **VWR**<sup>™</sup>

**VWR.COM**

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2021 Avantor, Inc. All rights reserved.

0221 LIT. NO. 080058W