



# MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System (and with all of the information required by the CPR)

Revision Date 07/05/2013

Version 1.2

## SECTION 1. Identification

### Product identifier

Product number 115934  
Product name Eosin B (bluish) (C.I.45400) for microscopy Certistain®

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses In vitro diagnostic reagent, Reagent for analysis

### Details of the supplier of the safety data sheet

Company EMD Millipore Corporation | 290 Concord Road, Billerica, MA 01821, United States of America | SDS Phone Support: +1-978-715-1335 | General Inquiries: +1-978-751-4321 | Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)

Emergency telephone 613-996-6666 CANUTEC (Canada)  
+1-703-527-3887 CHEMTREC (International)  
24 Hours/day; 7 Days/week

## SECTION 2. Hazards identification

### GHS-Labeling

Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008.

### Other hazards

None known.

## SECTION 3. Composition/information on ingredients

Chemical nature Dyestuff  
Formula  $C_{20}H_6Br_2N_2Na_2O_9$  (Hill)  
CAS-No. 548-24-3  
Molar mass 624.06 g/mol

Remarks WHMIS hazardous composition: No ingredients are hazardous according to the CPR criteria.

## SECTION 4. First aid measures

### Description of first-aid measures

#### Inhalation

After inhalation: fresh air.

## MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System (and with all of the information required by the CPR)

Product number

115934

Version 1.2

Product name

Eosin B (bluish) (C.I.45400) for microscopy Certistain®

---

### *Skin contact*

After skin contact: wash off with plenty of water. Remove contaminated clothing.

### *Eye contact*

After eye contact: rinse out with plenty of water.

### *Ingestion*

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

Never give anything by mouth to an unconscious person.

### **Most important symptoms and effects, both acute and delayed**

We have no description of any toxic symptoms.

### **Indication of any immediate medical attention and special treatment needed**

No information available.

---

## **SECTION 5. Fire-fighting measures**

### **Extinguishing media**

#### *Suitable extinguishing media*

Water, Carbon dioxide (CO<sub>2</sub>), Foam, Dry powder

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### **Special hazards arising from the substance or mixture**

Combustible.

Development of hazardous combustion gases or vapors possible in the event of fire.

Fire may cause evolution of:

nitrous gases, hydrogen bromide

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

#### *Further information*

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Avoid inhalation of dusts. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### **Environmental precautions**

Do not empty into drains.

### **Methods and materials for containment and cleaning up**

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

---

# MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System (and with all of the information required by the CPR)

Product number

115934

Version 1.2

Product name

Eosin B (bluish) (C.I.45400) for microscopy Certistain®

---

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

### Conditions for safe storage, including any incompatibilities

Tightly closed. Dry.

Store at +5°C to +30°C (+41°F to +86°F).

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

Contains no substances with occupational exposure limit values.

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### Hygiene measures

Change contaminated clothing. Wash hands after working with substance.

### Eye/face protection

Safety glasses

### Hand protection

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### Respiratory protection

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

## SECTION 9. Physical and chemical properties

Physical state

solid

Color

brown  
to  
green

Odor

weak

Odor Threshold

No information available.

---

## MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System (and with all of the information required by the CPR)

Product number	115934	Version 1.2
Product name	Eosin B (bluish) (C.I.45400) for microscopy Certistain®	

---

pH	ca. 6.2 at 10 g/l 77 °F (25 °C)
Melting point	No information available.
Boiling point	No information available.
Flash point	No information available.
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Vapor pressure	No information available.
Relative vapor density	No information available.
Relative density	No information available.
Water solubility	300 g/l at 77 °F (25 °C)
Partition coefficient: n-octanol/water	log Pow: -2.94 (calculated) (Lit.) Bioaccumulation is not expected (log Pow <1).
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Bulk density	ca.680 kg/m <sup>3</sup>

---

### SECTION 10. Stability and reactivity

#### Reactivity

The following applies in general to flammable organic substances and mixtures: in correspondingly fine distribution, when whirled up a dust explosion potential may generally be assumed.

#### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

# MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System (and with all of the information required by the CPR)

Product number

115934

Version 1.2

Product name

Eosin B (bluish) (C.I.45400) for microscopy Certistain®

---

## Possibility of hazardous reactions

Violent reactions possible with:

Strong oxidizing agents

## Conditions to avoid

no information available

## Incompatible materials

no information available

## Hazardous decomposition products

in the event of fire: See section 5.

---

## SECTION 11. Toxicological information

### Information on toxicological effects

#### *Likely route of exposure*

Eye contact, Skin contact, Ingestion

#### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

#### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

#### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

### Carcinogenicity

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

### Further information

Hazardous properties cannot be excluded but are unlikely when the product is handled appropriately.

Handle in accordance with good industrial hygiene and safety practice.

---

MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System (and with all of the information required by the CPR)

Product number

115934

Version 1.2

Product name

Eosin B (bluish) (C.I.45400) for microscopy Certistain®

---

**SECTION 12. Ecological information**

**Ecotoxicity**

No information available.

**Persistence and degradability**

No information available.

**Bioaccumulative potential**

*Partition coefficient: n-octanol/water*

log Pow: -2.94

(calculated)

(Lit.) Bioaccumulation is not expected (log Pow <1).

**Mobility in soil**

No information available.

---

**SECTION 13. Disposal considerations**

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

**SECTION 14. Transport information**

**Land transport (DOT)**

Not classified as dangerous in the meaning of transport regulations.

**Air transport (IATA)**

Not classified as dangerous in the meaning of transport regulations.

**Sea transport (IMDG)**

Not classified as dangerous in the meaning of transport regulations.

---

**SECTION 15. Regulatory information**

**United States of America**

**Canada**

**WHMIS Classification**

Not controlled under WHMIS (Canada).

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR.

---

## MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System (and with all of the information required by the CPR)

Product number	115934	Version 1.2
Product name	Eosin B (bluish) (C.I.45400) for microscopy Certistain®	

---

### Notification status

TSCA: All components of the product are listed in the TSCA-inventory.

DSL: All components of this product are on the Canadian DSL.

---

### SECTION 16. Other information

#### Training advice

Provide adequate information, instruction and training for operators.

#### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Revision Date 07/05/2013

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*