



Material	75164-300
Material description	L(+)-Potassium sodium tartrate tetrahydrate
Grade	ACS
CAS Number	6381-59-5
Molecular formula	C ₄ H ₄ KNaO ₆ ·4H ₂ O
Molecular mass	282,22

Characteristics

Specifications

Assay	99,0 --> 102,0 %
pH (25°C; 5 %)	6,0 --> 8,5
Heavy metals (as Pb)	Max. 5 ppm
Insoluble matter	Max. 0,005 %
Cl (Chloride)	Max. 0,001 %
NH ₄ (Ammonium)	Max. 0,002 %
PO ₄ (Phosphate)	Max. 0,002 %
SO ₄ (Sulphate)	Max. 0,005 %
Ca (Calcium)	Max. 0,005 %
Fe (Iron)	Max. 0,001 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada