



VWR ANALYTICAL

# Product Specification



Material BDH7540-4  
Material description Hydrogen peroxide 3% w/w

CAS Number 7722-84-1  
Molecular formula H<sub>2</sub>O<sub>2</sub>  
Molecular mass 34,01

Additional information

Characteristics	Specifications
Assay	3,00 --> 3,50 %
Appearance	Clear liquid
Stabilizer (Acetanilide)	Passes test
Prepared to formulation on file	Confirmed

Signature	Additional information
<p>We certify that this batch conforms to the specifications listed above.</p> <p>This document has been produced electronically and is valid without a signature.</p> <p>Alicia Bourdeau, Head of the Laboratory - Rouses Point VWR Chemicals, LLC, 3 Lincoln Blvd, Rouses Point, NY 12979-1045, USA</p>	<p>The weights and/or volumes used to prepare this product are NIST traceable. All balances used in the preparation of the product are calibrated weekly against NIST traceable weights. The balances are maintained and serviced on a regular basis by an outside certified company. All volumetric glassware used is NIST traceable and certified as meeting Class A specifications. Unless otherwise indicated above, all chemicals used in the preparation of this product are Reagent ACS grade.</p>