

VWR® Prefilled Histology Containers, 10% Neutral Buffered Formalin (NBF)

01. 10% FORMALIN
FIXATIVE FOR TISSUE
FIXATION

02. 20-120 ML
CONTAINERS ARE
95kPa-CERTIFIED

03. 180-1000 ML
CONTAINERS FEATURE
HEAT SEALS

04. ALSO AVAILABLE
IN BULK FORMULATION





VWR® PREFILLED FORMALIN HISTOLOGY CONTAINERS

- Each lot tested for formaldehyde concentration to ensure optimal fixation
- Manufactured in an ISO 13485 certified facility
- 20 mL through 120 mL sizes feature click-top and are 95kPa compliant
- 180 mL through 1000 mL sizes are heat-sealed to prevent leaking
- Strict process controls and lot tracking for full traceability
- Custom options available (minimum orders may apply)

Containers are half-filled with 10% neutral buffered formalin, for the collection and transport of tissue samples. Each container is made of high clarity virgin polypropylene and includes a polyethylene screw cap. Caps feature dual thread closures and a leak-resistant seal to protect specimens. All containers sized 20mL through 120mL are 95kPa compliant, for safe use in pneumatic tube systems.

Containers sized 20mL through 120mL feature a click-top closure. The first click assures users the seal is fully engaged. The second click ensures 95kPa compliance. This click-top design prevents leaking formalin containers and ensures safe specimen transport.

Containers sized 20mL through 120mL are packaged in transport packs with handles. This packaging design allows for easy distribution to other labs and doctors' offices. The inner packaging is made of chipboard material, and is safe for operating room use. Our large size prefilled containers (180mL through 1000 mL) are heat sealed to prevent leaking during shipping.

VWR® Bulk 10% Formalin is available in a 32 oz. bottle, gallon container, 2.5 and 5 gallon Cubitainer®, and a 55 gallon HDPE drum. Our cubitainers come complete with a spigot to allow easy filling of empty containers in the laboratory.

Size	Cat. No.	Unit
VWR Prefilled Formalin, Half Filled		
20 mL	16004-112	Cs. 128/
40 mL	16004-115	Cs. 48/
60 mL	16004-117	Cs. 54/
90 mL	16004-119	Cs. 48/
120 mL	16004-121	Cs. 48/
180 mL	16004-122	Cs. 50/
240 mL	16008-000	Cs. 24/
480 mL	16004-124	Cs. 24/
1000 mL	16004-130	Cs. 12/
VWR Bulk Formalin		
1 L	10790-714	Cs. 12/
1 gal.	89370-094	Cs. 4/
2.5 gal.	16004-126	Each
5 gal.	16004-128	Each
55 gal.	10790-716	Each



800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.