


FOR A LIMITED TIME

PURCHASE \$1000 of products and **GET UP TO \$500** additional products for **FREE!***

Pall's Microbiology QC products keep the 'Quality' in Quality Control



Improve efficiencies and flexibility with Pall's trusted microbiology QC solutions and expertise for water testing.

Pall offers a complete line up of Micro QC products including vacuum manifolds, portable dispensers & pumps, and a wide variety of membrane and media.

**Purchase a minimum of \$1000 of Pall Microbiology QC products in one order between January 1, 2022 to December 31, 2022 and get up to \$500 of FREE product. To redeem for your free product, visit vwr.com/promotions and search for promo code 5711. Fill in the online form before January 31, 2023. Be sure to indicate your choices of free product. Limit one redemption per customer. \$1000 purchase and \$500 credit must be for different products. Value of free products is based on VWR list price. Products Ordered or Free must be from list associated for this promotion. Promotion ends December 31, 2022.*

- Visit vwr.com/pall for information on qualifying products
- For further details contact your Avantor or Pall rep.

Ask me about Legionella, Cryptosporidia and Giardia testing!

FEATURED PRODUCTS


Sentino Pump



- Compact design – Small footprint frees valuable benchtop space
- Minimal Height – Convenient to use in a hood
- Disposable fluid path – Eliminate potential biofilm build-up

- Easy to clean – Wipe down with standard lab disinfectants

Magnetic Filter Funnel and Membranes



- No leak magnetic design allows one-handed operation
- Reusable and environmentally friendly
- Sturdy and safe. Polyphenylsulfone construction provides durability and added safety at less cost than most glass funnels

Top-Selling Products for Water Testing

Purchase \$1000 Worth of Products Below

Pall Laboratory Manifold

VWR Cat. No.	Pall Part No.	Description	Pkg
75981-548	4889	3 place Manifold	1/pkg
75981-554	4892	Standard Adapter	3/pkg
75981-552	4891	Sentino Adapter	3/pkg
75981-560	4959	Elongated Adapter	3/pkg
75981-550	4890	MicroFunnel Adapter	3/pkg
75981-556	4893	Coupling Device	1/pkg

Sentino Dispenser

VWR Cat. No.	Pall Part No.	Description	Pkg
10147-394	13184	Sentino Dispenser	1/pkg

Magnetic Filter Funnel

VWR Cat. No.	Pall Part No.	Description	Pkg
28143-550	4242	Magnetic Funnel 300 mL, 47 mm	1/pkg
28143-546	4247	Magnetic Funnel 150 mL, 47 mm	1/pkg

Sentino Microbiology Pump

VWR Cat. No.	Pall Part No.	Description	Pkg
10841-650	13186	Sentino Pump	1/pkg

Get up to \$500 Worth of the Following Products FREE

MicroFunnel Filter Funnels

VWR Cat. No.	Pall Part No.	Description	Pkg
28143-544	4800	GN6 Grid 0.45 µm, 100 mL	50/pkg
28143-560	4815	GN6 Grid 0.45 µm, 300 mL	20/pkg
28143-542	4803	Supor® Grid 0.2 µm, 100 mL	50/pkg
28143-556	4805	Metricel® Black Grid 0.45 µm, 47 mm, 100 mL	50/pkg
97058-848	FMFNL1050	Polycarbonate 0.4 µm, 47 mm, 100 mL	50/pkg
97058-850	FMFNL3020	Polycarbonate 0.4 µm, 47 mm, 300 mL	20/pkg

Media

VWR Cat. No.	Pall Part No.	Description	Pkg
28145-680	4302	Broth M-FC, 2 mL	50/pkg
28150-500	68105	Broth M-Endo, 2 mL	50/pkg

Forceps

VWR Cat. No.	Pall Part No.	Description	Pkg
30033-042	51147	Forceps	1/pkg

Membranes

VWR Cat. No.	Pall Part No.	Description	Pkg
28148-926	66278	GN6 Grid 0.45 µm, 47 mm	200/pkg
28147-979	66234	Supor Grid 0.2 µm, 47mm	200/pkg
28149-472	66585	Metricel Black 0.45 µm, 47mm	1000/pkg

Sentino Filter Funnels

VWR Cat. No.	Pall Part No.	Description	Pkg
10147-760	4870	Sentino Filter Funnel, 100 mL	100/pkg
10147-762	4871	Sentino Filter Funnel, 250 mL	80/pkg

Free product may also include any MicroFunnel or Sentino filter funnels and GN Metricel membranes with list prices below \$500



Ordering:
vwr.com

© 2022 Pall Corporation. Pall, Envirocheck, Metricel, MicroFunnel, Sentino, and Supor are trademarks of Pall Corporation. ® indicates a trademark registered in the USA.

2/22, GN20.0612ENVIRO



Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.