

Certificate of Quality



Part Number	MFLX95802-10
Lot Number	28574873
Part Description	TUBING SILICONE 3/16X7/16IN 25FT
Expiration Date	16DEC2026
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: FALTTR50-5
- The components/materials used in the production of the product were manufactured in a manner compliant with good manufacturing practices.
- This product is manufactured in the U.S.A.
- Meets or exceeds USP Class VI requirements.
- Meets or exceeds USP Class VI, USP 88 Class VI and USP 87 requirements, current revision.
- Meets or exceeds per 21 CFR 177.2600 "Rubber articles intended for repeated use may be safely used in producing, manufacturing, processing, preparing, treating, packaging, transporting, or holding food, subject to the provisions of this section."
- Also per 21 CFR 177.2600(G) "In accordance with good manufacturing practices finished rubber articles intended for repeated use in contact with food shall be thoroughly cleansed prior to their first use in contact with food."
- Listed materials are entirely non-toxic, non-hemolytic and non-pyrogenic.
- This material is not intended for invasive use.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.