

Corning® LSE™ Platform Rocker

Ideal for blotting and staining applications



The Corning LSE Platform Rocker is a compact yet high capacity rocker that is ideal for western blotting and general mixing applications. The adjustable rocking speed and fixed 7 degree tilt makes it optimal for staining fragile gels, as well as the more intensive action necessary for washing blots.

The speed of the Corning Platform Rocker can be adjusted from 2 to 60 rpm and the timer can be set for up to 2 hours or for continuous operation. The unit is safe to use in temperature controlled environments up to 70°C and is available in single or high capacity double platform configurations. Both of these rugged models come with a 5-year warranty.

Features and Benefits

- ▶ Low profile and compact size
- ▶ 11-pound load capacity
- ▶ Non-slip rubber pads to prevent sample containers from moving during operation
- ▶ Double platform model available for increased capacity



Specifications

Speed Settings	5 to 120 tilts/minute
Timer	0 to 2 hours or continuous
Motion	Rocking
Tilt Angle	Fixed 7°
Maximum Load	11 lbs. (5 kg)
Ambient Operation Range	+4°C to 70°C
Dimensions	12.8 x 14.2 x 5.7 in. (32.5 x 36 x 14.5 cm)
Weight	Single Platform: 11.7 lbs. (5.3 kg) Double Platform: 12.7 lbs. (5.8 kg)
Electrical	115V, 50/60Hz 230V, 50/60Hz
Warranty	5 years
Platform Size	11.8 x 11.8 in. (30 x 30 cm)

Ordering Information


Corning® LSE™ Platform Rockers

VWR Cat. No.	Corning Cat. No.	Description	Qty/Cs
76200-230	6703	Corning LSE Platform Rocker, single platform, 115V	1
76200-232	6704	Corning LSE Platform Rocker, single platform, 230V	1
76200-234	6705	Corning LSE Platform Rocker, double platform, 115V	1
76200-236	6706	Corning LSE Platform Rocker, double platform, 230V	1

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.

800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.