



RightCare™ NON-CONTACT INFRARED THERMOMETER

 Touch Free /
No Contact Required

 Measures
Temperature
within .5 Seconds

 Accurate /
Measures to +0.2°F

 Fever Alarm /
Automatic Alarm
Sound at 100.4°F



The RightCare™ Non-Contact Infrared Thermometer displays body temperature by measuring thermal radiation on the forehead. A built-in microprocessor and sensitive infrared sensor provide highly accurate and rapid results. No contact is required, making it a simple and safe way to take a temperature reading of everyone from a group of people to a sleeping baby. It can also measure the temperature of objects such as foods and liquids. Crisp LCD display can be seen in light or dark environments. Can store up to 32 readings. Includes 2 AA Batteries.

MPN: RCIRDТ

Each
/ Case of 50

CONTACT: Sales@ExpressMed.com | (877) 822-2719 Ext. 4343